

DONJOY®

GenuForce®

ELASTIC KNITTED KNEE SUPPORT
WITH KNEECAP PAD



enovis™

GenuForce®

Elastic knitted knee support with kneecap pad

GenuForce's® contoured knee pad helps protect and support the patella, while its flexible bilateral stays help provide relief and stabilization for the knee joint.

INTENDED USE/INDICATIONS:

- Support and compression following acute and chronic knee injury and mild knee osteoarthritis
- Proprioceptive and neuromuscular stimulation through compression
- Patellofemoral pain (knee pain)
- Mild knee instability
- Prevention of injury



ANATOMICALLY CONTOURED SILICONE KNEE PAD

- for additional stability of the patella, massaging effect, and to improve proprioception




FLEXIBLE BILATERAL STAYS

- help provide medial and lateral support



ORDERING INFORMATION

DESCRIPTION	PART NO.	Ⓕ THIGH CIRCUMFERENCE	Ⓖ CALF CIRCUMFERENCE	SIZE
	82-1750-1	37 - 40 cm 14½ - 15¾"	25 - 28 cm 9¾ - 11"	1 = XS
	82-1750-2	40 - 43 cm 15¾ - 17"	28 - 31 cm 11 - 11¾"	2 = S
	82-1750-3	43 - 46 cm 17 - 18"	31 - 34 cm 11¾ - 13½"	3 = M
GenuForce®	82-1750-4	46 - 49 cm 18 - 19¼"	34 - 37 cm 13½ - 14½"	4 = L
	82-1750-5	49 - 52 cm 19¼ - 20½"	37 - 40 cm 14½ - 15¾"	5 = XL
	82-1750-6	52 - 55 cm 20½ - 21¾"	40 - 43 cm 15¾ - 17"	6 = XXL
	82-1750-7	55 - 58 cm 21¾ - 22¾"	43 - 46 cm 17 - 18"	7 = XXXL
	82-1751-5	52 - 55 cm 20½ - 21¾"	37 - 40 cm 14½ - 15¾"	5+
GenuForce® Plus	82-1751-6	55 - 58 cm 21¾ - 22¾"	40 - 43 cm 15¾ - 17"	6+
	82-1751-7	58 - 62 cm 22¾ - 24½"	43 - 46 cm 17 - 18"	7+
	82-1751-8	62 - 65 cm 24½ - 25½"	46 - 49 cm 18 - 19¼"	8+



Ⓕ THIGH CIRCUMFERENCE
Measurement taken 15 cm above mid-patella

Ⓖ CALF CIRCUMFERENCE
Measurement taken 5 cm below mid-patella