



#### DONJOY®

# DonJoy ActyLight®

## Fit-and-forget ankle support

The convenience and comfort of the DonJoy ActyLight® ankle support means patients will forget they're even wearing it. Its step-in design and lockable lace system provide a fast, secure fit that's complemented by a pair of removable anatomical stays. Their simple adaptability enables both acute and chronic use, so ActyLight can support patients throughout their treatment.

### Intended use/Indications

DonJoy ActyLight is designed to provide support and stabilization to the ankle.

- · Initial management of chronic ankle instabilities
- · Mild to moderate sprains of the lateral ligament
- · Isolated syndesmotic injuries (Grade 1)
- · Return to sports, prevention of re-injury
- Ligament instability
- Supination prevention

## Features and benefits



Step-in design makes donning and doffing quick and easy



Secure lace-locking mechanism enables precise tightening, helps maintain a consistent fit, and allows quick release



Lightweight, breathable fabrics help the foot stay cool and dry to aid compliance



Figure-8 straps provide inversion/ eversion control



Anatomically contoured bilateral stays provide additional stability and protection and can be repositioned to accommodate the malleoli



Fits inside a trainer comfortably allowing patients to stay active



Modular design means the support can be adapted easily for acute or chronic use by removing the stays







ActyLight's step-in-sleeve design provides proprioceptive feedback and works with supportive elements to stabilize the ankle and restrict excessive range of movement.

- Bilateral semi-rigid supports provide mechanical stability to prevent ankle inversion/eversion
- Figure-8 straps provide further ankle inversion/eversion control
- Forefoot mesh made from breathable material to provide proprioceptive and compressive support for the foot



#### **ORDERING INFORMATION**

PART NO.	HEEL CIRCUMFERENCE	EU SHOE SIZE	SIZE
82-0054-1	≤ 27 cm	32 - 36	XS
82-0054-2	27 - 30 cm	36 - 39	S
82-0054-3	30 - 34 cm	39 - 43	М
82-0054-4	≥ 34 cm	43 - 47	L





Asia Pacific

1905, Tower 2, Grand Central Plaza

Shatin

HONG KONG

- 050 74054

E: marketing.ap@enovis.com

Canada

6485 Kennedy Rc

Mississauna

Ontario, L5T2W

CANADA

**T**: 1-866-866-503

E: canada.orders@enovis.com

Germany

Bötzinger Straße 90

79111 Freiburg

GERMANY

**T**: +49 (0) 180 1 676 333

E: kundenservice@enovis.com

South Africa

Unit 6, 2 on London, 2 London Circle

BrackenGate Business Park

Brackenfell, 7560

- -- (a) a: --- aa

E: info.southafrica@enovis.com

Australia – New Zealand

PO Box 6057

Frenchs Forest DC

**NSW 2086** 

AUSTRALIA

**r**: 1300 66 77 30

E: customerservice au@enovis com

...

Kleinhoefstraat 5, bus 39

Geel 2440

BELGIUN

**T**: +32 (0)14 248090

E: orders.bx@enovis.com

France

rue de Rethar

Centre Européen de Frê

64990 Mouguerre

FRANCE

**T**: +33 (0)5 59 52 86 90

E: sce.cial@enovis.com

Nordic

Aurmansgatan 126

SWEDEN

**T**: +46 40 39 40 00

E: info.nordic@enovis.com

Spain

Av.Cornella 144. 1-4

08950 Esplugues de llobrega

Barcelona

**T**: +34 93 480 32 02

E: ventas@enovis.com

UK & Ireland

1a Guildford Business Parl

Guildford, Surrey

002 0/

**T**: +44 (0)1483 459 659

E: ukorders@enovis.com

USA

2900 Lake Vista Drive

Lewisville

TX 75067

U.S.A.

**T**: 760 727 1280

E: customercare@enovis.com

#### enovis...

Creating Better Together\*

enovis-medtech.eu

1a Guildford Business Park, Guildford, Surrey, GU2 8XG, UNITED KINGDOM

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.