



enovis™

Creating Better Together™

CATALOGUE

enovis[™]

Creating Better Together[™]

CATALOGUE

WORKING WITH ENOVIS

Enovis® is a world leading orthopaedic company focused on helping people maintain and regain the joy of natural mobility. Offering a diverse portfolio of products and solutions to prevent injury, rehabilitate muscles, bones and joints and treat pain, we strive to help people live their lives to the fullest. We believe that our powerful brands, comprehensive range of products, extensive distribution network, unwavering focus on quality, innovation and customer service, and our strong relationships with orthopaedic and physical therapy professionals have all contributed to our leading market position. We have award winning manufacturing facilities across the globe with industry leading product design teams to continue the development and production of new products.

Our customer service team in South Africa is recognised in the industry as providing excellence in the support of our customers. Enovis Product Consultants throughout the South Africa provide detailed product and sales support.

ENOVIS CUSTOMER CARE

In addition to fast and accurate delivery, you can expect a friendly and knowledgeable Customer Service Department to take your order and answer any questions. Each of our Product Consultants serve and support an individual geographic territory, so you can count on expertise specific to your area. Should your area Product Consultant be unavailable, any of our other representatives will be equally helpful and able to serve your needs.

You can be assured your order will ship the same day, providing the item(s) requested are available and the order is placed before 3:00 PM. We look forward to serving you between the hours of 8.00^{AM}-4:30^{PM}.

HOW TO ORDER

To order from South Africa:

Landline: 021 276 2968

E-Mail: orders.zaf@enovis.com

QUESTIONS WE WILL ASK

Account Address: We will need to take the full delivery and invoice address.

Product: Please state the name, part number, description, quantities and size.

Practice name/number:

SPECIAL REQUIREMENTS

The majority of orders will ship within 24-72 Hours. Custom products and special orders (Replacement parts, special neoprene colours, etc.) are shipped within 6 weeks.

WARRANTY

All products have a limited warranty from the date of purchase against defects in materials and workmanship. Warranty varies with some products. Please contact Customer Services with specific questions.

RETURN POLICY

Before returning any products, please complete the Enovis Returns form and email it to our head office. A waybill will be generated and sent to you to be attached to your parcel with a copy of the invoice.

Please see page 106 for full terms and conditions.

OUR BRANDS

AIRCAST®

 chattanooga®

COMPEX 

DONJOY®




Dr.Comfort®

 exos®

PROCARE®

THE BEST FIT

We at Enovis have categorised our diverse portfolio into three areas, **SMART**, **STANDARD** and **VALUE**, to help you to identify the product best suited to your requirements.



THE CATEGORIES:

SMART Effective products that accelerate recovery with high compliance and simplified ease-of-use

STANDARD Everyday, essential products at a competitive price

VALUE Simple solutions for optimising the cost-benefit ratio under tight economic constraints

	SMART	STANDARD	VALUE
ACCELERATED RECOVERY	○	○	○
INNOVATIVE SOLUTIONS	○	○	○
PATIENT-CENTRIC DESIGN	○	○	
PATIENTS WITH ACTIVE LIFESTYLES	○	○	
CLINICALLY PROVEN	○	○	
RANGE OF INDICATIONS	○	○	○
TIME-PROVEN PRODUCTS	○	○	○
ECONOMY		○	○

CONTENTS

[THE BEST FIT]

THE BEST FIT.....	5
-------------------	---

[KNEE & LOWER LIMB]

4-POINTS-OF-LEVERAGE™	
BRACING SYSTEM.....	12
HOW DONJOY 4-POINTS-OF-LEVERAGE PROTECTS ACL.....	13
FOURCEPOINT HINGE TECHNOLOGY.....	13
KNEE BRACE LADDER	14
LIGAMENT DEFICIENCIES.....	15
MEASURING GUIDE	15
DEFIANCE® PRO ●	16
DEFIANCE COLOUR OPTIONS.....	17
DEFIANCE HINGE OPTIONS	17
FULLFORCE™ ●	18
ARMOR™ ●	18
STABILIZED KNEE SUPPORT ●	19
REDDIE® BRACE ●	19
QUICKFIT HINGED KNEE ●	19
PLAYMAKER XPERT ●	20
RIGID KNEE LIGAMENT BRACES	
ACCESSORIES	22
OSTEOARTHRITIS (OA)	23
ROAM OA™ ●	25
OA GO ●	25
OA DEFIANCE ●	26
OA ADJUSTER™ 3 ●	26
OA NANO™ ●	27
OA REACTION™ WEB ●	28
REACTION™ WEB ●	28
TRU-PULL® LITE ●	29
LATERAL J ●	29
INFRAPATELLAR BAND ●	29
POST OPERATIVE KNEE & HIP	30
X-ROM™ ●	31
X-ACT™ ROM LITE ●	32
PLAYMAKER X-ROM ●	32
COOL I-ROM ●	33
IMMO AT4 3V KNEE IMMOBILISER ●	33
ABDUCTION PILLOW ●	33
VERSAROM HIP™ ●	33

[FOOT]

ACTYTOE ●	36
REMEDY PRO™ ●	37
SQUARED TOE POST-OP SHOE ●	37
PODALUX II ●	37
PODAHEEL™ ●	38
MED/SURG SHOE ●	38
OFF-LOADING DIABETIC SHOE ●	39
AIRLIFT™ PTTD BRACE ●	39
SUPER-LITE A.F.O. ●	39
SILICONE HEEL CUPS ●	40
PODOUS BOOT ●	40
PLANTAR FASCIITIS NIGHT SPLINT ●	40
PODALIB® ●	40

[WALKER BOOTS]

AIRSELECT™ ELITE ●	44
AIRSELECT™ STANDARD ●	45
AIRSELECT™ SHORT ●	45
XP DIABETIC WALKER™ ●	46
XP DIABETIC WALKER REPLACEMENT PARTS ●	46
TUBE STRETCH SOCK ●	47
HEEL WEDGES ●	47
10° WEDGE KIT ●	47
MAXTRAX XCEL ROM ●	48
XCELTRAX™ TALL ●	48
XCELTRAX™ ANKLE ●	48
MINITRAX™ ●	49
SHOELIFT™ ●	49

[ANKLE]

ACTYLIGHT ●	52
ACTYFOOT™ ●	53
A60™ ANKLE SUPPORT ●	54
AIRSPORT+® ●	54
AIRHEEL™ ●	54
AIR-STIRRUP® UNIVERSE™ ●	54
AIR-STIRRUP® ●	55
STABILISING PRO ●	55
STABILISING SPEED PRO ●	55

[SPINE]

FORM FIT COLLAR ●	58
CALIFORNIA COLLAR ●	58
C1 COLLAR ●	59
ABDOSTRAP II ●	59
HERNIA BELT ●	59
PREMIUM PANEL ELASTIC BINDER ●	59
EXOS® FORM II ●	60
ISOFORM™ BACK BRACES ●	62
ISOFORM TLSO ●	64
ISOFORM POSTURAL EXTENSION ●	64
ACTISTRAP ●	65
COMFORTFORM BACK ●	65
IMMOSTRAP™ ●	65
MYBABYSTRAP ●	65
BODYSTRAP ●	66
SPINOSTRAP ●	66
OSTEOSTRAP ●	67

BEST FIT KEY

See page 5 for BEST FIT guide

- SMART
- STANDARD
- VALUE

CONTENTS

[SHOULDER & ARM]

SULLY® SHOULDER STABILISER ●	70
HUMERAL FRACTURE BRACE/OVER THE SHOULDER ●	70
HUMERAL STABILISING SYSTEM (HSS) ●	71
CROMAX ●	71
CLINICAL SHOULDER IMMOBILISER ●	71
HEALTHCARE BEAR ARM SLING ●	71
ULTRASLING QUADRANT™ ●	72
ULTRASLING® III ●	72
ULTRASLING PRO™ ●	73
QUICKFIT SHOULDER IMMOBILISER ●	73
X-ACT ROM ELBOW ●	74
PNEUMATIC ARMBAND ●	74
EPIFORCE STRONG ●	74

[WRIST & HAND]

WRIST BRACE WITH BOA® ●	78
RADIAL GUTTER FRACTURE BRACE WITH BOA® ●	78
HUMERAL FRACTURE BRACE WITH BOA® ●	78
SHORT ARM BRACE WITH BOA® OPEN THUMB ●	78
THUMB SPICA FRACTURE BRACE WITH BOA® ●	79
BOXER'S FRACTURE BRACE WITH BOA® ●	79
PAEDIATRIC SHORT ARM BRACE WITH BOA® - OPEN THUMB ●	79
SHORT THUMB SPICA II ●	79
EXTENDED SHORT THUMB SPICA ●	80
LONG THUMB SPICA II WITH BOA® ●	80
HAND BASED ULNAR GUTTER BRACE ●	80
HAND BASED RADIAL GUTTER BRACE ●	80
SHEET MATERIALS ●	81
BOA ADJUSTABLE TENSIONING KIT ●	81
STRAPS ●	81
TOP LOADING OVEN ●	81
RESPIFORM ●	82
RESPIFORM PLUS ●	82
DIGIFORM ●	82
DIGIFORM PLUS ●	82
CARPAFORM ●	84
POLYFORM WRIST ●	84
RHIZOFORM™ ●	85
ERGOFORM™ ●	85
COMFORTFORM WRIST ●	85
COMFORTFORM™ WRIST WITH ABDUCTED THUMB ●	85
QUICK-FIT WRIST ●	86
QUICK-FIT W.T.O ●	86
UNIVERSAL WRIST-O-PRENE ●	86
UNIVERSAL THUMB-O-PRENE™ ●	86
10" LACE UP ●	87
FROG SPLINTS / BASEBALL SPLINTS ●	87
STAXX PLASTIC FINGER SPLINTS ●	87

[HOT/COLD THERAPY]

CRYO/CUFF™ WRAPS & CUFFS ●	89
ICEMAN® CLEAR3™ ●	90
ICEMAN® WRAP-ON COLD PAD ●	90
CRYO/CUFF™ GRAVITY COOLER ●	91
FLEXIPAC ●	91

[COMPEX]

SP 8.0	94
SP 6.0	94
SP 4.0	94
MINI SET	94
FIXX 2.0 MASSAGER	95
AYRE	95

[VASCULAR]

VENAFLOW® ELITE ●	97
VENAPURE® ●	98
ANNIE BLACK ●	100
LANA BLACK ●	100
FLUTE BLACK ●	100
DIANE BLACK ●	100
MALORY BLACK ●	100
PROTECTOR BLACK ●	101
PERFORMANCE BLACK/GREY ●	101
PERFORMANCE BLACK ●	101
BRIAN BLACK ●	101
CLASSIC BLACK ●	101

[MISCELLANEOUS]

PRODUCT WARRANTY	103
TERMS AND CONDITIONS	104

BEST FIT KEY

See **page 5** for BEST FIT guide

- SMART
- STANDARD
- VALUE



DEFIANCE PRO

KNEE & LOWER LIMB

KNEE BRACING AT THE HEART OF ENOVIS HISTORY

While knee braces are a common sight, that hasn't always been the case. The modern knee brace, used to prevent injury, has only been around since 1967.

The late 60s were important for the evolution of knee braces. In 1967, a professor of physiology, Dr. Robert F. McDavid, invented the first lateral knee brace designed to prevent injury or re-injury. His brace provided lateral protection of the knee. Then during the 1967 season and Super Bowl III, Hall of Fame New York Jets' Quarterback Joe Namath played with a now-famous knee brace. Jack Castiglia of the Lenox Hill Brace Shop along with noted sports physician Dr. James Nicholas designed the brace. This allowed Namath to keep playing despite being plagued by knee problems. Namath's knee brace was so iconic that it is now on display at the Pro Football Hall of Fame.

During the 1970s, more prophylactic knee braces began being tested and tried by more NFL players in an effort to reduce career ending knee injuries. In 1978, Philadelphia Eagles' offensive line captain Mark Nordquist and lawyer Ken Reed founded a small company in a Carlsbad, California garage. They named the new company DonJoy® after their wives, Donna and Joy.

The first Enovis products were simple sleeves made of sewn-together neoprene that were pulled over the knee as well as the ankle and elbow joints for support. Fast forward to today where knee braces come in a variety of shapes and sizes, designed for all manner of knee injuries and prevention.

COLD THERAPY KNEE OPTIONS

For more details see pages 80-89



CRYO/CUFF™ KNEE



CRYO/CUFF™

4-POINTS-OF-LEVERAGE™ BRACING SYSTEM

The 4-Points-of-Leverage is a dynamic bracing system made up of a rigid cuff and strap configuration, whereby a posterior force is applied to the tibia, preventing anterior movement. This reduces the strain on the anterior cruciate ligament (ACL).

THE 4-POINTS-OF-LEVERAGE BRACE DESIGN:

- Significantly reduces ACL strain*
- Provides a constant dynamic load



POINT OF LEVERAGE #1

- The tibia is anchored by a cuff placed at the bottom of the posterior calf

POINT OF LEVERAGE #2

- The femur is anchored by a cuff placed at the top of the anterior thigh

POINT OF LEVERAGE #3


- A strap across the back of the lower thigh pushes the femur anteriorly

POINT OF LEVERAGE #4

- A strap placed on the anterior tibia applies an active constant load to prevent anterior tibia translation

*Fleming et al. (2000), The influence of functional Knee Bracing on the ACL Strain Biomechanics in Weight bearing and Nonweightbearing Knees, The American Journal of Sports Medicine, Vol. 6, p. 815-824

HOW DONJOY 4-POINTS-OF LEVERAGE PROTECTS ACL



Patients with an ACL deficient knee have an anterolateral instability that causes the knee to give way.

It is very difficult to control instability in an ACL deficient knee in a direct way using a knee brace.

Indirect stabilisation of the knee is achieved with the biomechanical 4 point design of the DonJoy brace.


The 4-point system is clinically proven to be the most effective way to prevent anterior translation of the tibia.

A DonJoy brace will provide a high level of stability for the ACL deficient knee by effectively reducing or preventing the anterior tibial translation through the use of the 4-points-of-leverage.

FOURCEPOINT HINGE TECHNOLOGY

Aggressive eccentric contractions of the quadriceps produce significant anterior tibial translation and can injure the ACL or permanently stretch out an ACL graft (DeMorat et al., AJSM, 2004). Additional studies show quadricep-intensive activities should be avoided at joint angles, particularly the last 30 degrees of extension, when strains are at a maximum. Enovis has developed a hinge to address maximising ACL and ACL graft protection in these "at risk" positions.

Small knee flexion angles at initial contact during landing tasks may contribute to ACL injuries. The FourcePoint Hinge can be used to control joint range of motion by increasing the knee flexion angle during landing, resulting in a decrease of anterior shear force on the ACL by 9-13% (Garret and Yu, AJSM, 2004).



The DonJoy FourcePoint utilises a spring mechanism in the hinge to apply gradually increasing resistance during knee extension. The resistance engages in the last 25 degrees of extension (relative to the extension stop).

The hinge has three levels of resistance with the option to turn off resistance if appropriate.

THE RESISTANCE SERVES THREE CRITICAL ROLES:

- First, it reduces the time spent near full extension or in the "at risk" position
- Second, it increases the posterior load applied by the brace on the tibia to prevent anterior tibial translation
- Third, it eliminates the extension shock felt when a patient extends into a 10 degree standard rigid stop


FourcePoint Hinge Technology works to enhance DonJoy's 4-Points-of-Leverage design by damping knee joint extension, which improves the mechanical performance of the brace and reduces shear forces at the knee. Addressing the "at-risk" position, the hinge resistance engages in the last 25 degrees of knee extension.

The FourcePoint hinge combined with the 4-Points-of-Leverage cuff and strapping design provides a more comfortable brace that reduces anterior shear forces at the knee. This is beneficial for patients as it reduces strain on the healing ACL graft, providing stability to an unstable, ACL deficient knee or for prophylactic use.

KNEE BRACE LADDER


Understand the product choices you have available to your knee injury patients from DonJoy.

ACL MCL LCL BRACING LADDER




*Playmaker Xpert for MCL and LCL instability only


CI BRACING LADDER



PCL BRACING LADDER



OA LADDER



LIGAMENT DEFICIENCIES

BRACE	LIGAMENT DEFICIENCIES SUPPORTED					BRACE OPTIONS/CHARACTERISTICS			
	ACL	MCL	PCL	LCL	COMBINED INSTABILITY (ACL & PCL)	HINGE	TYPE	STRAPPING OPTIONS	CALF CUFF OPTIONS
DEFIANCE PRO	○	○	○	○	○	Polycentric FourcePoint Drop-Lock OA	Rigid, customised	ACL, PCL or CI	Standard/Short
FULLFORCE	○					FourcePoint	Rigid	ACL	Standard
ARMOR	○	○	○	○	○	Polycentric	Rigid	CI	Short
QUICKFIT HINGED KNEE	○	○		○		Polycentric	Soft, sleeve	N/A	N/A
PLAYMAKER XPERT		○		○		Polycentric	Wraparound or sleeve	N/A	N/A
REDDIE BRACE		○		○		Polycentric	Soft, wraparound	N/A	N/A

MEASURING GUIDE

Measure at three points on the leg and use the chart to determine the size required.

Thigh 15cm above mid-patella (kneecap) · Knee, mid-patella (kneecap) · Calf 15cm below mid-patella (kneecap)

If one of the three measurements is a size up, select that size

For example, Thigh Size 2, Knee Size 3, Calf Size 2. Order Size 3

Note that for OA braces, the off-loading hinge is located on the opposite side of the affected compartment.

SIZE	THIGH CIRC. 15cm ABOVE MID-PATELLA	KNEE CENTRE CIRC. MID-PATELLA	CALF CIRC. 15cm BELOW MID-PATELLA
1 (XS)	33 – 39 cm	30.5 – 33cm	25.5 – 30.5cm
2 (S)	39 – 47 cm	33 – 35.5cm	30.5 – 35.5cm
3 (M)	47 – 53 cm	35.5cm – 38cm	35.5cm – 40.5cm
4 (L)	53 – 60 cm	38 – 43cm	40.5 – 47cm
5 (XL)	60 – 67 cm	43 – 48.25cm	47 – 50.75cm
6 (XXL)	67 – 75 cm	48.25cm – 53.25cm	50.75cm – 56cm
7 (XXXL)	75 – 81 cm	53.25 – 58.5cm	56 – 61cm

DONJOY

DEFIANCE[®] PRO



DonJoy's Defiance[®] has been the world's leading custom knee brace for over 20 years. Over 2 million Defiance braces have been prescribed to elite professional athletes, middle-school soccer players, and athletes of all shapes, sizes, ages and genders.

Now, Defiance PRO pushes our dominance even farther, boasting the greatest strength-to-weight ratio of any knee brace on the market.

Professional athletes choose Defiance PRO for the ultimate confidence and freedom to move.

INDICATIONS

Moderate to severe ACL, PCL, CI, MCL and LCL instabilities.

The OA Defiance is suitable for the active patient with osteoarthritis.



CHARACTERISTICS

- New low-profile design: The trusted high-performance carbon fiber frame now lighter and 5mm lower profile than Defiance Classic
- New C-6 liner material: Soft to the touch, moisture wicking and anti-microbial
- Less than 600g: Custom carbon fiber frame is built tough
- New silicone condyle pads: Comfortably help keep the brace in place
- Memory foam condyle pads: Contour to the knee (not shown)
- Swiveling internally mounted straps: Intimate fit around the anatomy for a lighter, more comfortable feel and improved anti-migration
- New soft touch strap ends: Flexible and strong hold keeps straps secure

ORDERING INFORMATION

SIZE	THIGH CIRC.	CALF CIRC.
Minimum	31.8 cm (12 ½")	26 cm (10 ¼")
Maximum	97 cm (38 ⅝")	66.5 cm (26 ¼")

Each Defiance is custom made. Measurements taken with CCMI Mark III

DEFIANCE PRO COLOUR OPTIONS



DEFIANCE HINGE OPTIONS



FOURCEPOINT ACL:

This patented hinge design helps protect the ACL by incorporating a dampening mechanism in the last 25° of extension to help keep the knee out of the at risk position and increase the posterior tibial load.

FOURCEPOINT PLUS:

The Plus hinge provides ACL protection benefits of the FourcePoint hinge with a reinforced steel outer plate ideally suited for the rigors of high-impact or collision sports

OA ADJUSTABLE HINGE:

The custom OA frame offers 1 to 7 degrees of off-load. The telescoping hinge helps provide additional off-loading up to 3 degrees. The memory foam condyle pads are included for increased patient comfort and compliance.

STANDARD LOW-PROFILE HINGE:

The polycentric hinge stabilises the knee while matching natural knee kinematics. Polycentric hinges also maximize the anterior/posterior and varus/valgus stabilising forces of the brace. This can be used for ACL, PCL, combined instability, MCL, or LCL indications.

DONJOY
FULLFORCE™



INDICATIONS

Moderate to severe ACL, MCL & LCL instabilities.
 For active lifestyles involving mid-range sports e.g. tennis, netball, cricket.

CHARACTERISTICS

- Enhanced 4-Points-of-Leverage with ACL strapping
- Available in 7 sizes, XS-XXXL
- Low profile FourcePoint hinge
- Strong, mouldable aircraft quality aluminium frame with carbon fibre look finish
- Internally mounted swivelling strap tabs to help accommodate leg movement
- Swooping thigh cuff for inner thigh clearance & bilateral use
- Optional short or standard calf length
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		DESCRIPTION
LEFT	RIGHT	
11-0265-x	11-0264-x	ACL, Standard Calf
13-8925		FourcePoint Allen Key

Sizing information is on page 15

DONJOY
ARMOR™



INDICATIONS

Moderate to severe ACL, PCL, CI, MCL and LCL instabilities. For active lifestyles involving extreme sports e.g. skiing.

CHARACTERISTICS

- Strapping configurations available for ACL, PCL or CI instabilities
- Available in 7 sizes, XS-XXXL
- Mouldable, medial swooping frame for custom fit and short calf for boot clearance
- Low profile hinges also offer flexion-extension control
- Strong aircraft quality aluminium frame to reduce brace weight and improve comfort
- Internally mounted strap tabs
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		DESCRIPTION
LEFT	RIGHT	
11-1030-x-26000	11-1029-x-26000	CI, Short Calf

Sizing information is on page 15



PROCARE
STABILIZED
KNEE SUPPORT



INDICATIONS

Ideal for mild patella tendonitis, chondromalacia or subluxations.

CHARACTERISTICS

- 3mm neoprene with medial/lateral spiral stays for added support
- Superior and inferior compression straps help prevent migration and provide additional support
- Universal neoprene buttress may be adjusted to allow for proper patella management and stabilisation

ORDERING INFORMATION

PART NUMBER	THIGH CIRC.	SIZE
79-82722	34cm-39cm	XS
79-82723	39cm-45.5cm	S
79-82725	45.5cm-52cm	M
79-82727	52cm-58.5cm	L
79-82728	58.5cm-64cm	XL
79-82729	64cm-71cm	XXL
79-82729-10	71cm-78.5cm	XXXL

Measurement taken 15cm above mid-patella



PROCARE
REDDIE® BRACE



INDICATIONS

Ideal for mild MCL and/or LCL sprains, mild arthritis, post arthroscopic meniscectomy, or mild global instabilities.

CHARACTERISTICS

- 3mm neoprene wraparound configuration suitable for the older patient
- Removable dual axis polycentric hinges provide medial/lateral support
- Popliteal cut-out for extra comfort in flexion
- Fits left or right knee
- Available in 8 sizes XS-XXXL
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	THIGH CIRC.	SIZE
79-82392	34cm-39cm	XS
79-82393	39cm-45cm	S
79-82395	45cm-52cm	M
79-82397	52cm-58cm	L
79-82398	58cm-64cm	XL
79-82399	64cm-71cm	XXL
79-82399-10	71cm-77cm	XXXL

Measurement taken 15cm above mid-patella

DONJOY
QUICKFIT HINGED KNEE



INDICATIONS

Protection, compression and support for knee sprains, mild to moderate ligament lesions, ligament instabilities, meniscus injuries, patella subluxation, laxity or misalignment, mild arthritic pain, return to mobility or sports activities.

CHARACTERISTICS

- Patella donut maintain patella in correct position.
- Donut can be trimmed to provide tailored patella misalignment correction.
- Removable aluminium polycentric hinges offer medial / lateral stability.
- Encased hinges for protection from lateral impacts (ideal for sports use).
- Soft, anti-pinch posterior lycra panel enhances comfort.
- Breathable, form-fitting 5 mm spacer mesh material promotes airflow.
- Fits Left or Right knee.

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
82-0010-1	26,6cm-33cm	XXS
82-0010-2	33cm-39cm	XS
82-0010-3	39cm-47cm	S
82-0010-4	47cm-53cm	M
82-0010-5	53cm-60cm	L
82-0010-6	60cm-67cm	XL
82-0010-7	67cm-75cm	XXL
82-0010-8	75cm-81cm	XXXL

Measurement taken 15cm above mid-patella



Playmaker Expert
Wraparound

DONJOY

PLAYMAKER XPERT



INDICATIONS

For mild medial/lateral support of the knee. For use in warmer environments.



CHARACTERISTICS

- Advanced high tenacity nylon/spandex material promotes a cooling effect while delivering a lighter and slimmer profile
- Removable polycentric hinges and aluminium uprights
- Popliteal cut-out for extra comfort in flexion
- Fits left or right knee
- 6 month warranty
- Latex-free

ORDERING INFORMATION

PART NUMBER		THIGH CIRCUMFERENCE	SIZE
WRAPAROUND	PULL-ON SLEEVE		
82-0023-1	82-0022-1	33 - 39cm	XS
82-0023-2	82-0022-2	39 - 47cm	S
82-0023-3	82-0022-3	47 - 53cm	M
82-0023-4	82-0022-4	53 - 60cm	L
82-0023-5	82-0022-5	60 - 67cm	XL
82-0023-6	82-0022-6	67 - 75cm	XXL
82-0023-7	82-0022-7	75 - 81cm	XXXL

Measurement taken 15cm above mid-patella



RIGID KNEE LIGAMENT BRACES ACCESSORIES



FULL EXTENSION STOP

PART NO.	SIZE
11-0398-9-00000	U



LYCRA UNDERGARMENT

PART NO.	SIZE
11-0075-x-00000	S-XXL



SPORTS COVER

PART NO.		SIZE
BRACE LENGTH - STANDARD	BRACE LENGTH - SHORT	
11-0016-x-06000	11-0103-x-06000	S-XXL



SUSPENSION STRAP KIT

PART NO.	SIZE
11-0543-1-06000	XS
11-0543-2-06000	S-M
11-0543-3-06000	L-XXL



NEOPRENE UNDERGARMENT

PART NO.	SIZE
11-0122-x	XS-XXL



AIR CONDYLE KIT

PART NO.	SIZE
25-1923-0-06000	U



IMPACT GUARD

PART NO.	SIZE
11-1006-x-06000	S-XXL



CONDYLE KIT

PART NO.	SIZE
25-0271-0-00000	U



FOAM CONDYLE KIT

PART NO.	SIZE
11-1107-0-00000	U



PINCH GUARD

PART NO.		SIZE
ARMOR BRACE	DEFIANCE BRACE	
11-1095-9-06000	11-1367-9-06000	U


OSTEOARTHRITIS (OA)



We are at a unique stage in our healthcare system where patients are being diagnosed with knee OA at a younger age, and there is increasing demand for conservative treatment options. OA is the fourth leading disability in the world with 8.75 million people in the UK seeking treatment for OA*. This combined with a growing active aging population sees the healthcare system not being able to continue to offer total knee replacements as a primary intervention.

OA knee bracing is a conservative management alternative.

*Reference-Arthritis Research UK



DONJOY®

THE NEXT GENERATION OF
OFFLOADER BRACES FROM DONJOY





NEW



INDICATIONS

For patients with symptoms, physical findings or radiographic evidence of moderate to severe unicompartmental knee osteoarthritis, post-operative meniscal tears only where unicompartmental unloading is indicated, or varus/valgus knee conditions that may benefit from unicompartmental unloading.

CHARACTERISTICS

- Easy mag-align clips.
- Set and forget technology
- Boa** dial for just the right fit
- Slim profile hinge
- Anti-migration technology
- Patent-pending comfort straps
- Hand Washable
- Backed by DonJoy®
- Easy to wear
- Only 6 off-the-shelf SKUs
- Condyle pad is included for increased patient comfort and compliance
- Suited for patients with active lifestyles or athletes.

ORDERING INFORMATION

PART NUMBER		THIGH CIRC.	KNEE CENTER	CALF CIRC.	SIZE
RIGHT MEDIAL/ LEFT LATERAL	LEFT MEDIAL/ RIGHT LATERAL				
11-1726-1-06000	11-1725-1-06000	33 - 47 cm	30.5 - 35.5 cm	25.5 - 35.5 cm	XS/S
11-1726-3-06000	11-1725-3-06000	47 - 59.5 cm	35.5 - 35.5 cm	35.5 - 45.5 cm	M/L
11-1726-5-06000	11-1725-5-06000	59.5 - 81 cm	43.2 - 58.4 cm	45.5 - 61 cm	XL/XXL/ XXXL

Note: Hinge is on the SAME side of affected compartment



NEW



INDICATIONS

DonJoy OA GO is designed to relieve pain, improve proprioception and knee function by supporting the body in guiding and stabilising the knee joint. For the use of mild to moderate single compartment osteoarthritis as well as varus-valgus knee conditions that may benefit from unicompartmental unloading.

CHARACTERISTICS

- Breathable, stretchable and antibacterial fabric enable long periods of use.
- Rear Lycra section to prevent bunching
- Hybrid sleeve/wrap design for custom fit
- Easy-grip pull tabs for easy don and doff of the brace
- Dual individual BOA fitting systems offer quick and easy custom relief when fastening
- Anatomically contoured condyle pads
- Anti-slip bands to prevent migration
- 4 flexible shells to provide optimal pressure distribution and anatomical fit

ORDERING INFORMATION

PART NUMBER		THIGH CIRCUMFERENCE	SIZE
RIGHT MEDIAL/ LEFT LATERAL	LEFT MEDIAL/RIGHT LATERAL		
82-7428-1	82-7429-1	33 – 39 cm	XS
82-7428-2	82-7429-2	39 – 47 cm	S
82-7428-3	82-7429-3	47 – 53 cm	M
82-7428-4	82-7429-4	53 – 60 cm	L
82-7428-5	82-7429-5	60 – 67 cm	XL
82-7428-6	82-7429-6	67 – 75 cm	XXL

Hinge is on opposite side of affected compartment
 Measurement taken 15 cm above mid-patella
 Available on request

DONJOY

OA DEFIANCE



INDICATIONS

Unicompartmental OA with or without combined ACL, PCL, Combined, MCL and LCL instabilities and activities of daily living.

CHARACTERISTICS

- Customised design, with a wide range of colour choices
- Reinforced carbon composite frame
- Durable support for moderate to severe levels of osteoarthritis.
- This polyaxial hinge comes with a telescoping condyle pad
- Allows patient-controlled additional unloading on the base load of 1 to 7° built into the frame.
- Condyle pad is included for increased patient comfort and compliance
- Suited for patients with active lifestyles or athletes.

ORDERING INFORMATION

SIZE	THIGH CIRCUMFERENCE	CALF CIRCUMFERENCE
Minimum	31.8 cm (12 1/2")	26 cm (10 1/4")
Maximum	97 cm (38 1/8")	66.5 cm (26 1/4")

Each Defiance is custom made. Measurements taken with CCMI Mark III

DONJOY

OA ADJUSTER™ 3



INDICATIONS

For patients with moderate to severe unicompartmental OA, needing significant off-load and ligament stability in a tough durable frame.

CHARACTERISTICS

- Nano MAG™ technology on the off-loading side, combined with the durability of the aluminum frame, provides lightweight function with optimal strength
- Clinically Proven-DonJoy's Adjuster technology provides adjustable off-loading with the highest level of condylar separation
- 4-Points-of-Leverage system-provides clinically proven ligament stability combined with low-profile hinges and double upright design allow a contoured, personalised fit for all levels of activity
- Available in 7 sizes, XS-XXXL
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

PART NUMBER		AFFECTED COMPARTMENT
LEFT	RIGHT	
11-1591-x	11-1590-x	Medial Compartment
11-1593-x	11-1592-x	Lateral Compartment

Note: Adjustment hinges are on opposite side of affected compartment. Sizing information is on page 15

PART NUMBER	DESCRIPTION
13-3120-0-06000	Replacement Key
25-1376-0-00000	Stop kit, extension 20°, 30°, 45° flexion 60°, 75° and full lockout



DONJOY
OA NANO™

Combining Nano MAG™, a new hi-tech material, with DonJoy's proven adjuster technology, OA Nano is the world's lightest functional OA knee brace.

INDICATIONS

For mild to moderate unicompartmental OA, the OA Nano is ideal for patients who want to stay active and move freely.



CHARACTERISTICS

- The flexibility of the magnesium frame providing an intimate fit and increased comfort with Nano MAG (on the loading side) provides strength for off-loading
- DonJoy's clinically proven three point off-loading technology shifts the compressive loads from the degenerative area to the healthy compartment

- 4 Points of Leverage and contoured frame design enhances brace suspension for all levels of activity
- Available in 7 sizes, XS-XXXL
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

PART NUMBER		AFFECTED COMPARTMENT
LEFT	RIGHT	
11-1215-x	11-1214-x	Medial Compartment
11-1217-x	11-1216-x	Lateral Compartment

Note: Adjustment hinges are on opposite side of affected compartment. Sizing information on page 15

PART NUMBER	DESCRIPTION
13-3120-0-06000	Replacement Key
25-1376-0-00000	Stop kit, extension 20°, 30°, 45° flexion 60°, 75° and full lockout

DONJOY

OA REACTION™ WEB



INDICATIONS

Suitable for early to mild unicompartmental osteoarthritis of the knee & patella. The ideal "first-time" choice, allowing patients to get accustomed to wearing a brace and making it easier to transition into a more rigid OA brace when the time comes. OA Reaction Web is also a good option for patients wishing to use it as a lighter alternative to their usual rigid brace.

CHARACTERISTICS

- Four way stretch wrap-around design
- Silicone web technology
- Low profile, lightweight, discreet
- Single aluminium upright
- Silicone condyle pad
- Fitted undersleeve
- Available in 7 sizes, XS-XXXL
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	AFFECTED COMPARTMENT
11-7426-x	Right Medial/Left Lateral
11-7427-x	Left Medial/Right Lateral

Note: Off-loading hinge is on opposite side of affected compartment. Sizing information is on page 15

DONJOY

REACTION™ WEB



INDICATIONS

Chondromalacia patella, Osgood-Schlatter disease, patella OA, mild bicompartmental OA, quadriceps and patella tendonitis/tendinosis.

CHARACTERISTICS

- Flexible web design absorbs shocks and peak loads by dispersing energy across the knee
- The web design dynamically stabilises the patella on all sides with dual axis hinges for additional support
- Fitted under-sleeve for comfort
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE	SIZE
11-0215-2	33cm-46cm	XS/S
11-0215-3	47cm-59cm	M/L
11-0215-4	60cm-75cm	XL/XXL
11-0215-5	>75cm	XXXL

Measurement taken 15cm above mid-patella

ACCESSORIES	DESCRIPTION	SIZE
11-0541-x	Replacement Undersleeve	S to XXXL

Sizing information is on page 15



INDICATIONS

Recommended for acute mild to moderate anterior knee pain, support for patella instability, lateral subluxation and maltracking.

CHARACTERISTICS

- Patented dual durometer (silicone) buttress system provides consistent corrective force on the patella
- Bifurcated strap dynamically pulls the patella into proper alignment
- Semi-flexible removable nylon hinges included
- 21cm long support giving maximum power to pull patella into alignment
- Fabric is hypoallergenic, breathable, and lightweight
- 6 month warranty

ORDERING INFORMATION

PART NUMBER		THIGH CIRCUMFERENCE	SIZE
LEFT	RIGHT		
11-0261-1	11-0260-1	33cm-39cm	XS
11-0261-2	11-0260-2	39cm-47cm	S
11-0261-3	11-0260-3	47cm-53cm	M
11-0261-4	11-0260-4	53cm-60cm	L
11-0261-5	11-0260-5	60cm-67cm	XL
11-0261-6	11-0260-6	67cm-75cm	XXL
11-0261-7	11-0260-7	75cm-81cm	XXXL

Measurement taken 15cm above mid-patella



DONJOY
LATERAL J



INDICATIONS

The tubular Lateral "J" buttress is ideal for supporting lateral patellar subluxation, dislocations or lateral patellofemoral malalignment.

CHARACTERISTICS

- "J"-shaped buttress for patellar stabilisation
- Lateral to medial strapping for additional stabilisation
- Drytex material for patients with allergies or irritations to neoprene

ORDERING INFORMATION

PART NUMBER		THIGH CIRC.	SIZE
LEFT	RIGHT		
81-06602	81-06592	33 - 39cm	XS
81-06603	81-06593	39 - 47cm	S
81-06605	81-06595	47 - 53cm	M
81-06607	81-06597	53 - 60cm	L
81-06608	81-06598	60 - 67cm	XL
81-06609	81-06599	67 - 75cm	XXL

Measurement taken 15cm above mid-patella



AIRCAST
INFRAPATELLAR
BAND



INDICATIONS

Osgood Schlatter disease; patella tendonitis.

CHARACTERISTICS

- An aircell focuses compression on the patella tendon, helping to decrease stress at the tibial tubercle
- Fits left or right knee
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	CIRCUMFERENCE	SIZE
08A-B	25.5cm-43cm	U

Measurement taken just below the kneecap

POST OPERATIVE
KNEE & HIP





DONJOY
X-ROM™



INDICATIONS

To aid immobilisation and protected range of motion associated with ACL, PCL, LCL, and MCL surgeries and meniscal repairs



CHARACTERISTICS

- Quickly fit patients right out of the box
- Precisely control ROM with settings in 10° increments from -10° to 90° extension and -10° to 120° flexion
- Quick-lock design: Easily lock at 10° increments from -10° to 30°
- Strap, tab and buckle design make post-op application and adjustment fast and simple
- Easily accommodate patients from 125cm - 198 cm in height

- Four independently telescoping sliders with easy push button control settings help enhance brace placement, stability and comfort while avoiding the incision and peroneal nerve
- Bendable hinge bars help enhance varus/valgus contouring
- Secure contact points
- Strong, lightweight design helps reduce brace migration
- Anti-microbial pads

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
-------------	-------------	------

11-2181-9	X-ROM™ Post-Op Knee Brace	U*
-----------	---------------------------	----

*Length adjustability: 48 to 68.3cm

DONJOY

X-ACT™ ROM LITE



INDICATIONS

Designed for immobilisation and protected range of motion associated with post-traumatic treatment, ACL, PCL, LCL and MCL surgeries, patella realignment and meniscal or articular repairs.

CHARACTERISTICS

- Economy version of the popular X-Act ROM knee brace with easy to use strap system
- Adjustable uprights can be contoured to the patient's leg to provide varus and valgus support
- Extension settings from -10° to 90°
- Flexion settings from -10° to 120°
- Easy quick-lock at 0°, 15° and 30°
- Full circumferential foam
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
-------------	-------------	------

11-2161-9	X-Act ROM Lite	U*
-----------	----------------	----

*Length adjustability: 46cm to 65cm *Maximum thigh circumference: 76cm

DONJOY

PLAYMAKER X-ROM



INDICATIONS

For immobilisation and protected range of motion associated with ACL, PCL, LCL, and MCL surgeries and meniscal repairs. Can be used for young children and adolescences.

CHARACTERISTICS

- Precise ROM Protection. Control ROM with settings in 10-degree increments from -10 degrees to 90 degrees extension and -10 degrees to 120 degrees flexion.
- Quick-lock Design. Easily lock at 10-degree increments from -10 to 30 degrees.
- Comfortable Fit. High-tenacity (Hi10™) nylon/spandex fabric offers unmatched comfort and fit.
- Breathable spacer mesh fabric. Lightweight spacer fabric on front and back panels promotes a "cool" environment support.
- Fast, easy application. Moulded strap tabs allow quick and easy application and strap adjustments to deliver a secure and comfortable fit.
- Anti-microbial Coating. Neutralises odour while increasing comfort and durability.
- Customisable design. Trimmable, adjustable strap and buckle design allows off-the-shelf and customisable applications.

ORDERING INFORMATION

PART NUMBER	THIGH CIRC.	SIZE
-------------	-------------	------

11-0021-1	30.5cm-33cm	XS
-----------	-------------	----

11-0021-2	33cm-35.5cm	S
-----------	-------------	---

11-0021-3	35.5-38cm	M
-----------	-----------	---

11-0021-4	38cm-43cm	L
-----------	-----------	---

11-0021-5	43cm-48.5cm	XL
-----------	-------------	----

11-0021-6	48.5cm-53.3cm	XXL
-----------	---------------	-----

11-0021-7	53.3cm-58.5cm	XXXL
-----------	---------------	------



DONJOY
COOL I-ROM



INDICATIONS

For warmer environments, the Cool IROM provides immobilisation and precise range of motion control for ACL, PCL, MCL and LCL surgeries.

CHARACTERISTICS

- 10° incremental range of motion control
- Mouldable hinge bars allow for varus/valgus positioning
- Hinge bar breakdown capability
- 4 straps and abbreviated foam at key control points
- Integrated strap/hinge bar design
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	LEG INSEAM
11-0175-2-06000	Cool IROM Short	<80cm
11-0175-4-06000	Cool IROM Long	>80cm



DONJOY
IMMO AT4 3V KNEE
IMMOBILISER



INDICATIONS

Pre or post-operative total immobilisation of the knee joint.

CHARACTERISTICS

- Lightweight (less than 600 grams) universal 3 panel immobiliser
- Malleable & removable bilateral aluminium stays and a posterior stay
- Perforated foam for enhanced comfort
- Internal anti-migration strap
- Fits left or right leg
- Contains Latex

ORDERING INFORMATION

PART NUMBER	LENGTH	SIZE
82-43-30	30cm (12")	U/Paediatric
82-43-01	40cm (16")	U
82-43-02	50cm (20")	U
82-43-03	60cm (24")	U
82-43-04	65cm (26")	U

Measure inseam of leg minus 20cm (8")



PROCARE
ABDUCTION PILLOW



INDICATIONS

Designed for use following hip surgery when immobilisation or post-operative positioning is required.

CHARACTERISTICS

- All foam construction for maximum support
- Soft contact closure foam laminate straps

ORDERING INFORMATION

PART NUMBER	LENGTH x WIDTH	SIZE
79-90173	45cm x 30cm	S
79-90175	58cm x 38cm	M
79-90177	66cm x 46cm	L



DONJOY
VERSAROM HIP™



INDICATIONS

DonJoy® VersaROM™ Hip is indicated for use with post-op patients who have undergone a surgical procedure and need an orthosis until the component is integrated and soft tissue heals.

CHARACTERISTICS

- Integrated belt compression with cinching straps
- Malleable hip and thigh frames
- Precise ROM Adjustments
- Range of Motion allows controlled hip flexion and extension adjustments in 15° increments
- Telescoping Height Adjustment
- Laterally placed abduction paddle with three-point pressure
- Customised Fit
- Lightweight, slim-profile design

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
11-3335	DonJoy® VersaROM™ Hip

FOOT



ACTYTOE

FOOT BRACING



FOOT

THE BEST FIT GUIDE

POST-OP SHOES		
SMART	STANDARD	VALUE
Podalux II	Remedy Pro	Med/Surg Shoe
-	Squared Toe Post-Op Shoe	-
-	Off-Loading Diabetic Shoe	-
-	Podaheel	-

SPLINTS		
SMART	STANDARD	VALUE
AirLift PTTD Brace	PodaLib	Plantar Fasciitis Night Splint
-	Actytoe	Superlite A.F.O.
-	Podous Boot	Silicone Heel Cup



AIRCAST
ACTYTOE



INDICATIONS

- Pain relief in patients with symptomatic Hallux Valgus (bunion)
- Mild or moderate mal-positioning of the big toe
- Maintaining position after Hallux Valgus (bunion) surgery



CHARACTERISTICS

- Low-profile design
- Locking clip for splint immobilization (0° and 30° dorsi- and plantar flexion)
- Memory foam cushion
- Adjustable correction angle (0° and 10°)
- Adjustable, thin and soft straps
- Anatomically contoured splint
- Breathable and lightweight materials

- For day and night wear
- Suitable for walking
- Available in 3 sizes (S, M, L)
- For left or right foot
- Can be worn with wide footwear
- Not made with natural rubber latex

ORDERING INFORMATION

PART NUMBER	UK SHOE		SIZE
	MENS	WOMENS	
82-04S	< 4	< 3.5	S
82-04M	3.5 - 8	3 - 7.5	M
82-04L	7.5 - 11.5	7 - 11	L



INDICATIONS

Reducing forefoot pressure, forefoot trauma, post-surgical healing and diabetic ulcers or wounds.

CHARACTERISTICS

- 15° wedge design reduces forefoot pressure by shifting weight to midfoot and heel
- Easy forefoot closure
- Durable outer sole, padded liner and adjustable ankle strap
- Seamless, padded heel
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE		SIZE
	MENS	WOMENS	
79-81722	<4	<4.5	XS
79-81723	4 - 6	4.5 - 6.5	S
79-81725	6 - 8	6.5 - 8.5	M
79-81727	8 - 10	8.5 - 10.5	L
79-81728	>10	>10.5	XL



PROCARE
SQUARED TOE
POST-OP SHOE



INDICATIONS

Post-operative treatment of hallux valgus, toe deformities and protocols where maximum immobilisation and protection of the forefoot are required.

CHARACTERISTICS

- Squared design allows additional room for bulky dressings
- High ankle strapping for ease of ambulation and pressure reduction
- Rigid rocker sole with skid resistant sole
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE		SIZE
	MENS	WOMENS	
79-81232	2 - 3.5	3 - 4.5	XS
79-81233	4 - 5.5	5 - 6.5	S
79-81235	6 - 7.5	7 - 8.5	M
79-81237	8 - 10.5	9 - 11.5	L
79-81238	>10.5	>11.5	XL



DONJOY
PODALUX II



INDICATIONS

Post-op, post-trauma foot and toe applications and for patients requiring extra stability, more specifically after hallux valgus and ensuring the foot, toes and forefoot are held in a proper position during rehabilitation.

CHARACTERISTICS

- Wide, solid skid resistant sole for extra stability
- Removable rigid foot plate for additional firmness and arch support
- Lined with breathable micro-mesh textile for patient comfort
- Two contact closure straps for easy application and adjustment
- Toe piece included to cover up the toes if necessary
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	SHOE MENS	SHOE WOMENS	SIZE
82-0042-1	2.5 - 3.5	3 - 4	S
82-0042-2	4 - 5	4.5 - 5.5	M
82-0042-3	6 - 6.5	6 - 7	L
82-0042-4	7 - 7.5	7.5 - 8	XL
82-0042-5	8 - 10.5	8.5 - 10.5	XXL
82-0042-6	11 - 13	11 - 13	XXXL

DONJOY

PODAHEEL™



INDICATIONS

Reducing pressure on the heel and rearfoot. Treatment of lesions, pressure sores, diabetic ulcers or wounds and for post-surgical healing following bone spurs, fractures or stress fractures.

CHARACTERISTICS

- 10° wedge design places the foot in plantar flexion shifting weight to the forefoot
- Wide forefoot closure strap accommodates dressings
- Durable outer sole, and soft padded liner
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE MENS	UK SHOE WOMENS	SIZE
82-0031-1	4 - 6	3.5 - 5.5	S
82-0031-2	6 - 8	5.5 - 7.5	M
82-0031-3	8 - 11	7.5 - 10.5	L
82-0031-4	>11	>10.5	XL

Available on request

PROCARE

MED/SURG SHOE



INDICATIONS

Post-op, post-trauma applications where some mobility of the foot is allowed.

CHARACTERISTICS

- Loop/lock contact closure straps
- Comfortable contoured heel with foam mesh upper and semi-flexible skid resistant sole
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE MENS	SIZE
79-81132	<4	XS
79-81133	4 - 8	S
79-81135	8 - 11	M
79-81137	11 - 14	L
79-81138	>14	XL

PART NUMBER	UK SHOE SIZE WOMENS	SIZE
79-81142	<1	XS
79-81143	1 - 3.5	S
79-81145	3.5 - 5.5	M
79-81147	5.5 - 7.5	L

PROCARE

OFF-LOADING DIABETIC SHOE



INDICATIONS

Pressure relief during the treatment of diabetic ulcers.

CHARACTERISTICS

- Tri-laminate IMPAX grid insole designed to redistribute load away from ulcers with easy to remove IMPAX foam blocks
- IMPAX grid allows for a continuous plantar surface
- Lightweight breathable material enhances patient's comfort and promotes compliance
- Squared toe box and wide footbed to leave ample room for wound dressings
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE MENS	UK SHOE SIZE WOMENS	SIZE
79-81513	3.5 - 5.5	4 - 6	S
79-81515	6 - 8	6.5 - 8.5	M
79-81517	8.5 - 10.5	9 - 10	L
79-81518	11 - 13.5	>10.5	XL

ORDERING INFORMATION - REPLACEMENT INSOLE

PART NUMBER	UK SHOE SIZE MENS	UK SHOE SIZE WOMENS	SIZE
79-81523	3.5 - 5.5	4 - 6	S
79-81525	6 - 8	6.5 - 8.5	M
79-81527	8.5 - 10.5	9 - 10	L
79-81528	11 - 13.5	>10.5	XL



DONJOY

AIRLIFT™ PTTD BRACE



INDICATIONS

- Phase I and initial phase II of posterior tibial tendon dysfunction, first symptoms of flat foot.

CHARACTERISTICS

- Adjustable arch aircell for individual support and comfort
- Anatomically designed shells for secure ankle stability
- Integrated inflation pump
- Step-in design and simple two strap application
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER		UK SHOE MENS	UK SHOE WOMENS	SIZE
LEFT	RIGHT			
02PSL-1	02PSR-1	<6	<6.5	S
02PML-1	02PMR-1	6-10	6.5-10.5	M
02PLL-1	02PLR-1	>10.5	>11	L



PROCARE

SUPER-LITE A.F.O.



INDICATIONS

Correction of drop foot.

CHARACTERISTICS

- Low profile
- Trim to fit ankle orthosis
- Female (narrower) and male (wider) sizes
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
79-81150	Right, Male	U
79-81151	Left, Male	U
79-81160	Right, Female	U
79-81161	Left, Female	U



DONJOY
SILICONE HEEL
CUPS



INDICATIONS

Plantar fasciitis, heel spurs and Achilles tendonitis.

CHARACTERISTICS

- Moulded medical-grade soft silicone to reduce shock and vibration loads to joints
- Softer centre in heel region to reduce strain to area sensitive to pain
- Sold in pairs

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
81-81103	4.5 - 8.5	3 - 7	S/M
81-81107	>8.5	>7.5	L/XL



PROCARE
PODOUS BOOT



INDICATIONS

Prevent lower extremity disorders associated with trauma or immobility including: pressure necrosis, ankle contractures and post-operative hip abduction or rotation.

CHARACTERISTICS

- Lightweight, quality plastic shell with breathable poly/pile liner
- Adjustable toe extension helps protect toes from bedding
- Anti-rotation bar may be repositioned to assist with internal or external rotation
- Latex free
- Replacement liners available as a pack of 3

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
79-90550	<9.5	<9	U
79-90551	>9.5	>9	XL
REPLACEMENT LINERS		PACK SIZE	SIZE
79-90560		Pack of 3	U
79-90561		Pack of 3	XL



PROCARE
PLANTAR FASCIITIS
NIGHT SPLINT



INDICATIONS

Controlled application of a prolonged stretch on the plantar fascia and Achilles tendon.

CHARACTERISTICS

- Lightweight, sturdy one-piece flexible shell with a terry liner
- Adjustable bilateral dorsiflexion straps
- Padded calf and ankle straps for additional comfort
- Includes optional toe wedge for added stretch of the plantar fascia
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MENS	WOMENS	
79-97753	<6	<5	S
79-97755	6 - 8.5	5 - 8	M
79-97757	8.5 - 11	8 - 10.5	L
79-97758	11 - 13.5	>10.5	XL



DONJOY
PODALIB®



INDICATIONS

Drop foot, neurological weakness or paralysis of the foot-lifting muscles, helps prevent lateral imbalance without restricting movement, improving ambulation.

CHARACTERISTICS

- Quick closure buckles helps to adjust and secure the tension of the elastic strap and foot flexion angle
- Adjustable elastic strap with metal hook adapts to all shoe types with laces, buckles or eyelets
- Optional foot cuff for barefoot use or swimming
- Anatomical ankle cuff with two aluminium stays for improved support
- Posterior foam pad for comfort on the Achilles Tendon area
- Contains latex

ORDERING INFORMATION

PART NUMBER	ANKLE CIRC.	ANKLE HEIGHT	SIZE
82-0013-1	13cm - 19cm	10cm	XS
82-0013-2	13cm - 20cm	13cm	S
82-0013-3	20cm - 27cm	13cm	M
82-0013-4	27cm - 34cm	13cm	L





AIRSELECT ELITE

WALKER BOOTS



WALKER BOOTS

THE BEST FIT GUIDE

	WALKERS		
	SMART	STANDARD	VALUE
WALKERS	AirSelect Elite	AirSelect Standard/Short	XcelTrax Tall/Ankle
	Minitrax	Maxtrax Xcel ROM	-
	XP Diabetic walker	-	-

COLD THERAPY ANKLE OPTIONS

For more details see pages 80-89



CRYO/CUFF® ANKLE



AIRCAST

AIRSELECT™ ELITE



AirSelect puts comfort, control and healing in the hands of the patient with the most advanced pneumatic walking boot available. Beneath its sleek appearance are patented technological advances that work together to enhance treatment and improve outcomes. The integrated inflation system lets patients inflate each aircell individually, to maximise comfort and minimise oedema.

INDICATIONS

Specifically designed for stable fractures of the lower leg, foot, and ankle; severe ankle sprains; and post-operative use. Clinically proven effective in healing fractures and managing oedema.



CHARACTERISTICS

- The most advanced walking boot, engineered to provide the ultimate in protection, comfort, and oedema control
- Overlapping posterior and distal duplex aircells provide intermittent compression, reducing oedema 3 times faster
- Pre-inflated front panel aircell and custom inflation system provide a full-contact fit
- Enclosed toe cover provides protection from hazards and the elements
- Latex free
- Fits left or right leg
- 6 month warranty

PART NUMBER	SHOE SIZE	
	UK	SIZE
01EP-XS	<3	XS
01EP-S	3 - 6	S
01EP-M	6 - 9	M
01EP-L	9 - 12	L
01EP-XL	>12	XL

AIRCAST

AIRSELECT™ STANDARD



INDICATIONS

Stable fracture of foot/ankle/lower leg, severe ankle sprain, post-operative immobilisation.

CHARACTERISTICS

- Semi-rigid shell that supports the limb while providing protection
- Two adjustable distal aircells that provide compression and support
- Integrated inflation system for ease of use
- Lab-tested rocker sole encourages natural gait
- Latex free
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER	SHOE SIZE	SIZE
	UK	
01EF-XS	<3	XS
01EF-S	3 - 6	S
01EF-M	6 - 9	M
01EF-L	9 - 12	L
01EF-XL	>12	XL

AIRCAST

AIRSELECT™ SHORT



INDICATIONS

Metatarsal fractures, stress and stable fracture of the foot, bunionectomy, hallux valgus, soft tissue injury, severe ankle sprain.

CHARACTERISTICS

- The Short walking boot has a durable, semi-rigid shell that helps support the limb while providing protection
- The pre-inflated front panel and two customisable distal aircells provide compression and support
- Integrated inflation system for ease of use
- SoftStrike technology absorbs and dissipates shock, while the lightweight and lab-tested rocker sole encourages natural gait
- Latex free
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER	SHOE SIZE	SIZE
	UK	
01ES-XS	<3	XS
01ES-S	3 - 6	S
01ES-M	6 - 9	M
01ES-L	9 - 12	L
01ES-XL	>12	XL

AIRCAST

XP DIABETIC WALKER™



INDICATIONS

Stable foot and/or ankle fracture; severe ankle sprain; post-operative; Charcot foot; ulcer in malleolus, hind-foot and heel.

CHARACTERISTICS

- Specifically designed to meet the needs of the higher risk diabetic patient with a high rocker sole for maximum off-loading of the diabetic foot
- Dual density insole to help eliminate pressure points
- Graduated compression to help reduce oedema
- Regulated aircell inflation with the specialised hand bulb with pressure gauge
- Latex free
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SHOE SIZE	SIZE
	UK	
01PD-S	4 - 5.5	S
01PD-M	6 - 8	M
01PD-L	>8	L
01PD-XL	>10.5	XL

AIRCAST

XP DIABETIC WALKER REPLACEMENT PARTS



INDICATIONS

Replacement parts for use with XP Diabetic Walker System.

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
01PG	Replacement Pump with Pressure Gauge	U
01R-S	Diabetic Replacement Kit (includes 1 sock and 1 insole)	S
01R-M	Diabetic Replacement Kit (includes 1 sock and 1 insole)	M
01R-L	Diabetic Replacement Kit (includes 1 sock and 1 insole)	L
01R-XL	Diabetic Replacement Kit (includes 1 sock and 1 insole)	XL



AIRCAST
TUBE STRETCH
SOCK



INDICATIONS

To be used with Aircast walker boots to create a barrier between the boot and the skin, to absorb any sweat and decrease friction.

CHARACTERISTICS

- Cotton sock for use with AirSelect Elite, Short and Standard
- Fits left or right foot
- Latex free

ORDERING INFORMATION

PART NUMBER	SIZE
0129A	U



AIRCAST
HEEL WEDGES



INDICATIONS

Designed to be used with the Aircast walkers. Indicated for post-achilles tendon surgery, the heel supports can be placed on top of each other to achieve the desired heel rise.

CHARACTERISTICS

- Supplied as a pack of five heel wedges of varying size to suit most walker boots
- Up to three wedges can be combined in a single walker to achieve the following levels of plantarflexion:
 - - 3 heel wedges: 22 degrees
 - - 2 heel wedges: 16 degrees
 - - 1 heel wedge: 10 degrees
- Latex free

ORDERING INFORMATION

LEFT	RIGHT	WALKING BRACE SIZE
01K-L	01K-R	S/M/L
01K-XLL	01K-XLR	XL

AIRCAST
10° WEDGE KIT



INDICATIONS

The Aircast 10° wedges are designed: to comfortably position the ankle joint in plantar flexion during achilles tendon rehabilitation for use with the aircast airselect standard and elite walker boots

CHARACTERISTICS

- Customized height adjustments - the kit contains 3 wedges:
 - 3 wedges = angle of 30° of plantar flexion
 - 2 wedges = angle of 20° of plantar flexion
 - 1 wedge = angle of 10° of plantar flexion
- Lightweight but strong hollow design to keep additional weight to a minimum
- Double-sided adhesive tape included with each wedge for easy application and secure fixation
- Bilateral: fits either foot
- Latex free
- Sustainable design made from recycled material

ORDERING INFORMATION

PART NUMBER	AIRSELECT WALKER BOOT SIZE
82-0110	XS-S-M
82-0210	L-XL

MAXTRAX XCEL ROM



INDICATIONS

Promotes early range of motion for achilles tendon repair for controlled rehab and stable fractures of the lower leg.

CHARACTERISTICS

- Adjustable ROM hinge: allows ROM settings in 7.5° increments between 45° plantar-flexion and 30° dorsi-flexion with locking option
- Unique pre-shaped double upright design to take pressure of the injury at both the malleoli and lower leg
- Edgeless rocker sole is specially engineered to help promote a natural, stable gait
- Shock-absorbing sole reduces the impact of heel-strike during ambulation, minimising further irritation to the injury site
- Padded removable liners for cool comfort during rehabilitation and easy cleaning and replacement
- Adjustable hook-and-loop straps

ORDERING INFORMATION

MAXTRAX XCEL ROM	UK SHOE SIZE	SIZE
TALL		
82-104-XS	<3.5	XS
82-104-S	4 - 6.5	S
82-104-M	7 - 10	M
82-104-L	10 - 12	L
82-104-XL	>12	XL



PROCARE XCELTRAX™ TALL



INDICATIONS

Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

CHARACTERISTICS

- Lightweight design
- Unique pre-shaped ergonomic frame helps to modify the level of pressure on the injury at both the malleoli and lower leg
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE	SIZE
79-95492-1	<3.5	XS
79-95493-1	4 - 6.5	S
79-95495-1	7 - 10	M
79-95497-1	10 - 12	L
79-95498-1	>12	XL

AVAILABLE ON REQUEST



PROCARE XCELTRAX™ ANKLE



INDICATIONS

Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

CHARACTERISTICS

- Short version
- Lightweight design
- Unique pre-shaped ergonomic frame helps to modify the level of pressure on the injury at both the malleoli and lower leg
- Fits left or right leg

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE	SIZE
79-95502-1	<3.5	XS
79-95503-1	4 - 6.5	S
79-95505-1	7 - 10	M
79-95507-1	10 - 12	L
79-95508-1	>12	XL

AVAILABLE ON REQUEST

PROCARE
MINITRAX™



INDICATIONS

Ideal support and comfort for your smallest patients following acute ankle sprains, metatarsal fractures, stress and stable fractures of the lower leg.

CHARACTERISTICS

- Child friendly design elements to make wearing a boot less intimidating
- Lightweight design features anatomically shaped uprights to facilitate fit
- Low profile rocker sole features a durable, non-slip tread
- Easy to apply and adjust with a simple strapping and liner system
- Fits left or right leg

ORDERING INFORMATION

PART	AGE	FOOT	LEG	SHOE	SIZE
79-95510	1 - 2 ½	<15cm	<23cm	<7.5	S
79-95520	2 ½ - 5	15-17.5cm	<28cm	7.5 - 10	M
79-95530	5 - 8	17.5-20cm	<36cm	10.5 - 1 adult	L



SMALL

MEDIUM

PROCARE
SHOELIFT™



INDICATIONS

Simple, fast and easy to apply, the Shoelift™ uses an innovative design that helps to equalize a patient's limb length and reduce body strain while walking in a cast or walking boot.

CHARACTERISTICS

- Easily removable
- Recommended for securely fastened closed toe shoe or sneaker
- Available in three different sizes
- Fits left or right foot

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE		SIZE
	MALE	FEMALE	
79-95603	7-8.5	6.5 - 9.5	S
79-95605	9-11.5	9.5-12.5	M
79-95607	12+	12.5+	L



ANKLE

ACTYFOOT

ANKLE BRACING

For more than 30 years, medical professionals worldwide have trusted Aircast® for technological breakthroughs in patient care. Rooted in scientific methods, each Aircast product is developed using the concept of "functional management." Aircast set a new standard of care for sprains and other injuries with patented technology and the use of graduated pneumatic compression.

THE BEST FIT GUIDE

PREVENTIVE (CHRONIC ANKLE INSTABILITY)		
SMART	STANDARD	VALUE
-	A60	Stabilising pro
-	Actylight	Stabilising speed pro
-	ActyFoot	-

ANKLE SPRAIN GRADE I		
SMART	STANDARD	VALUE
-	Actylight	-

ANKLE SPRAIN GRADE II		
SMART	STANDARD	VALUE
AirSport+	ActyFoot	Air-Stirrup Universe
Air-Stirrup	-	-

ANKLE SPRAIN GRADE III		
SMART	STANDARD	VALUE
AirSelect Elite	AirSelect Standard	-

COLD THERAPY ANKLE OPTIONS

For more details see pages 80-89



CRYO/CUFF® ANKLE

NEW



DONJOY ACTYLIGHT



INDICATIONS

DonJoy ActyLight is designed to provide support and stabilisation to the ankle.



CHARACTERISTICS

- Step-in design
- Secure lace-locking mechanism enables precise tightening and allows quick release
- Lightweight, breathable fabrics help the foot stay cool and dry to aid compliance
- Figure-8 straps provide inversion/eversion control
- Anatomically contoured bilateral stays provide additional stability and protection and can be repositioned to accommodate the malleoli
- Fits inside a trainer comfortably
- Modular design means the support can be adapted easily for acute or chronic use by removing the stays

PART NUMBER	HEEL CIRC	UK SHOE SIZE	SIZE
82-0054-1	≤ 27 cm	<3	XS
82-0054-2	27 - 30 cm	3-6.5	S
82-0054-3	30 - 34 cm	6.5-9	M
82-0054-4	≥ 34 cm	9-12	L



AIRCAST
ACTYFOOT™



INDICATIONS

Designed to provide support, stabilisation and partial immobilisation of the ankle. Acute lateral ligament injuries. Acute ankle sprains and ligament ruptures. Prevention of further injuries in pre-operative chronic ankle instability. Protection of ankle joint



CHARACTERISTICS

- Three integrated straps provide proprioceptive feedback and work together with semi-rigid elements to stabilise the ankle and restrict excessive range of movement.
- Detachable lateral stay allows brace to be modified for use in both acute and chronic injury phases
- Low profile. Fits in most shoes or trainers and can be worn over a sock for additional comfort
- Can be fully opened to allow simple, secure application
- Calf strap provides additional stability to the ankle joint to help improve overall balance and control
- Cross strap further helps minimise the risk of ankle inversion
- Forefoot strap allows wearer to adjust compression as required to help manage swelling

ORDERING INFORMATION

PART NUMBER		HEEL CIRC.	UK SHOE SIZE	SIZE
LEFT	RIGHT			
82-02A-SL	82-02A-SR	<30cm	>5	S
82-02A-ML	82-02A-MR	30cm-34cm	5-8	M
82-02A-LL	82-02A-LR	>34cm	8-10	L
82-02A-PL	82-02A-PR	≥ 36 cm	>10	XL

Measurement taken across heel (calcaneus) to front of foot (talus)



AIRCAS
A60™ ANKLE
SUPPORT



INDICATIONS

Chronic instability, prophylaxis.

CHARACTERISTICS

- A 60° stabiliser guards against rollover and helps prevent sprains
- Breath-O-Prene® fabric keeps wearer cool and dry
- Sleek, anatomic design for a comfortable bulk-free fit
- Simple single strap application and adjustment
- Latex free
- Now also available in white
- 6 month warranty

ORDERING INFORMATION

PART NUMBER		UK SHOE SIZE	SIZE
LEFT	RIGHT		
02TSL	02TSR	<6	S
02TML	02TMR	6 - 8	M
02TLL	02TLR	>8	L



AIRCAS
AIRSPORT+®



INDICATIONS

Acute ankle sprain, post-operative (tendon rupture and ankle fracture), chronic instability, prophylaxis.

CHARACTERISTICS

- Boa® lacing system provides adjustable, uniform compression without creating pressure points
- Foam-filled aircell cushions which support the ankle
- Anterior cross strap and forefoot wrap provide additional compression and stability to the joint
- Contains latex

ORDERING INFORMATION

PART NUMBER		UK SHOE SIZE	SIZE
LEFT	RIGHT		
82-0266-1	82-0265-1	<3	XS
82-0266-2	82-0265-2	3 - 5	S
82-0266-3	82-0265-3	6 - 8	M
82-0266-4	82-0265-4	9 - 11	L
82-0266-5	82-0265-5	>11	XL



AIRCAS
AIRHEEL™



INDICATIONS

Chronic and acute Achilles tendonitis, calcaneal pain, plantar fasciitis.

CHARACTERISTICS

- Pulsating pneumatic compression minimises swelling and enhances circulation
- Lightweight materials and low profile design increase comfort
- Simplified design speeds application
- Fits left or right foot
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	UK SHOE SIZE	SIZE
09AS	<6	S
09AM	6 - 8	M
09AL	>8	L



AIRCAS
AIR-STIRRUP®
UNIVERSE™



INDICATIONS

Ankle sprains and chronic ankle instability.

CHARACTERISTICS

- Patented Duplex aircell system to reduce swelling and pain by enhanced circulation
- Universal semi-rigid shell design fits left or right ankle
- Universal size
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SIZE
02E	U



AIRCAST
AIR-STIRRUP®



INDICATIONS

Acute ankle distortion grade I, II & III, post-operative application, chronic ankle instability.

CHARACTERISTICS

- Semi-rigid, anatomically designed shells for protection, comfort and the prevention of inversion/eversion
- Patented Duplex aircell system to enhance circulation and reduce swelling
- Pre-inflated aircells for easy application
- Streamlined to fit in shoes for early protected weight-bearing
- Latex free
- 6 month warranty

ORDERING INFORMATION

PART NUMBER		PATIENT HEIGHT	BRACE LENGTH	SIZE
LEFT	RIGHT			
02JL	02JR	Small children	15cm	XS (Paed.)
02CL	02CR	<157cm (<5'2")	22cm	S
02BL	02BR	157cm-167cm (5'2"-5'6")	23cm	M
02AL	02AR	>162cm (>5'4")	27cm	L



DONJOY
STABILISING PRO



INDICATIONS

Chronic ankle sprains and prophylactic use.

CHARACTERISTICS

- Improved low profile design reduces the risk of ankle injury by up to 61%
- Wide (7.5cm) circumferential elastic strap for compression and stability
- Removable medial/lateral stays for additional support and stability
- Fits left or right ankle
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	HEEL/ANKLE CIRC.	SIZE
11-3234-0-06000	23cm - 25cm	XXS
11-3234-1-06000	25cm - 28cm	XS
11-3234-2-06000	28cm - 30cm	S
11-3234-3-06000	30cm - 33cm	M
11-3234-4-06000	33cm - 36cm	L
11-3234-5-06000	36cm - 38cm	XL
11-3234-6-06000	38cm - 41cm	XXL
11-3234-7-06000	41cm - 43cm	XXXL

Measurement taken from base of the heel around the ankle whilst wearing a sports sock

DONJOY
STABILISING
SPEED PRO



INDICATIONS

Chronic ankle sprains and prophylactic use.

CHARACTERISTICS

- Improved low profile design reduces the risk of ankle injury by up to 61%
- Wide (7.5cm) circumferential elastic strap for compression and stability
- Removable medial/lateral stays for additional support and stability
- Speed lacing system for fast application
- Fits left or right ankle
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	HEEL/ANKLE CIRC.	SIZE
11-3235-0-06000	23cm - 25cm	XXS
11-3235-1-06000	25cm - 28cm	XS
11-3235-2-06000	28cm - 30cm	S
11-3235-3-06000	30cm - 33cm	M
11-3235-4-06000	33cm - 36cm	L
11-3235-5-06000	36cm - 38cm	XL
11-3235-6-06000	38cm - 41cm	XXL
11-3235-7-06000	41cm - 43cm	XXXL

Measurement taken from base of the heel around the ankle whilst wearing a sports sock

SPINE



SPINOSTRAP

SPINE BRACING



THE BEST FIT GUIDE

SPINAL SUPPORT		
SMART	STANDARD	VALUE
Exos Form back braces	IsoFORM Postural Extension	Premium Panel Elastic Binders
IsoFORM TLSO	Actistrap	Comfortform back
Osteostrap	Immostrap	Hernia belt
My Babystrap	IsoFORM SIO	Bodystrap
Abdostrap	IsoFORM LO/LO+	Spinostrap
-	IsoFORM LSO/LSO+	-

NECK		
SMART	STANDARD	VALUE
-	C1 Collar	Formfit collar
-	-	California collar

SPINE

COLD THERAPY BACK OPTIONS

For more details see pages 80-89



CRYO/CUFF™ BACK/HIP/RIB

FORM FIT COLLAR



INDICATIONS

This collar provides semi-rigid support for use in soft tissue damage.

CHARACTERISTICS

- Contoured for anatomical fitting
- Comfortable and secure fit

ORDERING INFORMATION

PART NUMBER	HEIGHT X LENGTH	NECK CIRC.	SIZE
79-83043	7,5 x 47 cm	11 – 16 cm	S
79-83045	10 x 51 cm	13 – 18 cm	M
79-83047	11,5 x 57 cm	15 – 20 cm	L

CALIFORNIA COLLAR



INDICATIONS

Ideal for cervical spine immobilisation and extraction.

CHARACTERISTICS

- Pre-formed, lightweight, two-piece polyethylene foam with rigid plastic supports front and rear
- Collar maintains head in neutral position
- May be used to provide support in combination with a cervical or scoliosis cast

ORDERING INFORMATION

PART NUMBER	HEIGHT	NECK CIRC.	SIZE
79-83123**	5,7 cm	25,5 - 33 cm	S
79-83125**	5,7 cm	33 - 40,5 cm	M
79-83127**	5,7 cm	40,5 - 48 cm	L
79-83133*	8,3 cm	25,5 - 33 cm	S
79-83135*	8,3 cm	33 - 40,5 cm	M
79-83137*	8,3 cm	40,5 - 48 cm	L
79-83143**	10,8 cm	25,5 - 33 cm	S
79-83145**	10,8 cm	33 - 40,5 cm	M
79-83147**	10,8 cm	40,5 - 48 cm	L
79-83153**	13,3 cm	25,5 - 33 cm	S
79-83155**	13,3 cm	33 - 40,5 cm	M
79-83157**	13,3 cm	40,5 - 48 cm	L

Height: sternal notch to chin *Items in Stock **Available on request

NEW

DONJOY C1 COLLAR



INDICATIONS

Designed for cervical sprains and strains, post injury rehabilitation or post rigid collar use.

CHARACTERISTICS

- Design supports the cervical spine in a neutral position
- Low density foam collar
- Anatomically preformed at the chin and ear level for better fit
- 100% cotton covering with contact closure strap
- Available in navy blue

ORDERING INFORMATION

PART NUMBER	NECK CIRC.	SIZE
HEIGHT = 7.5CM		
CC10PPAB	21cm - 25cm	Paediatric
CC10P0AB	25cm - 28cm	XXS
CC10P1AB	29cm - 32cm	XS
CC10P2AB	33cm - 37cm	S
CC10P3AB	38cm - 42cm	M
CC10P4AB	43cm - 47cm	L
CC10P5AB	48cm - 51cm	XL



DONJOY HERNIA BELT



INDICATIONS

Provides adjustable, gradual compressive support to help reduce inguinal hernias and alleviate discomfort.

CHARACTERISTICS

- Easy to adjust, soft waist band and perineal straps provide comfortable compression
- Removable foam compression pads for gentle, gradual compression
- Treats right, left or double hernias
- Comfortable, lightweight nylon material
- Latex free

ORDERING INFORMATION

PART NUMBER	WAIST CIRC.	SIZE
81-97903	70 cm - 80 cm	S
81-97905	80 cm - 95 cm	M
81-97907	95 cm - 110 cm	L
81-97908	110 cm - 120 cm	XL



DONJOY ABDOSTRAP II



INDICATIONS

Abdominal support for temporary abdominal weakness, post-partum, post-op use following abdominal surgery.

CHARACTERISTICS

- Moderate compressive elastic material
- 2 Flexible abdominal stays for increased abdominal support
- Anterior panel: trimmable textile for probe or catheter insertion
- 2 Flexible dorsal stays
- Removable additional abdominal strap included for enhanced support

ORDERING INFORMATION

PART NUMBER			WAIST CIRC.	SIZE
18 CM*	25 CM	32 CM*		
82-0256-1	82-0257-1	82-0258-1	56 cm - 68 cm	S
82-0256-2	82-0257-2	82-0258-2	68 cm - 82 cm	M
82-0256-3	82-0257-3	82-0258-3	82 cm - 98 cm	L
82-0256-4	82-0257-4	82-0258-4	98 cm - 116 cm	XL
82-0256-5	82-0257-5	82-0258-5	116 cm - 136 cm	XXL

Available on request for the 32cm



PROCARE PREMIUM PANEL ELASTIC BINDER



INDICATIONS

Compression and support for post natal or post abdominal surgery, abdominal strains and weakness.

CHARACTERISTICS

- Constructed with panel elastic and soft cotton type lining for patient comfort
- Universal style support with contact closure
- Multi panel design helps prevent panel rolling
- Available in 22.5cm (9") 3 panel and 30cm (12") 4 panel heights

ORDERING INFORMATION

PART NUMBER		WAIST CIRC.	SIZE
9" 3-PANEL	12" 4-PANEL		
79-89070	79-89090	76-114 cm (30"-45")	S
79-89071	79-89091	114-157 cm (45"-62")	M
79-89210	79-89220	157-188 cm (62"-74")	L
79-89270	79-89280	188-213 cm (74"-82")	XL
79-89250	79-89260	213-254 cm (82"-100")	XXL

EXOS
EXOS® FORM II



V-STX™ MATERIAL

A patented, thermoformable material used with the posterior, anterior and lordotic inserts, allows for true customization and ideal patient fit.



COMFORTCORE™ FOAM

Conforms to the unique contours of a patient's body, adjusts instantly to individual movements.



DIAMOND2 GRID™

Uniquely designed construct, delivers patients' a more durable, longer-wearing brace.



BOA® FIT SYSTEM

The micro-adjustable BOA® Fit System allows for a fast, effortless, precision fit through a customisable and seamless level of support.



Exos Form II is designed with patented thermoformable VST-X™ material, Comfort CORE™ Foam that contours to the patient's body, Diamond Grid construct for a more durable and longer-wearing brace and Boa® Technology's advanced closure system.


The end result is a spine brace that delivers superior decompression and support, a quick and easily customizable fit, and provides patient comfort and satisfaction by reshaping form, fit and function.

FEATURES AND BENEFITS

- With a simple, lightweight, low profile design, the Exos Form™ II range is designed with patient comfort in mind without sacrificing functionality
- Adjustable (angle) belt wings accommodate varying patient body structures for optimal fit
- The modular step-up/step-down design enables you to adjust and configure your brace as you progress through rehabilitation and recovery
- A semi-universal design minimizes inventory management requirements

COMMON INDICATIONS

SACRO	T9-S1	T9-S1
621	631	637



Sacroiliac Pain	•		
Acute and Chronic Lower Back Pain		•	•
Lower Back Sprains / Strains		•	•
Lumbar Disc Displacement		•	•
Osteoporosis		•	•
Disc Herniation and Degeneration		•	•
Post-Operative Laminectomy		•	•
Post-Operative Disectomy		•	•
Spondylolisthesis		•	•
Post-Operative Fusion		•	•
Spinal Stenosis		•	•
Stable / Non-Displaced Fractures		•	•



EXOS® FORM II 637

- Provides a higher degree of spinal support and relief for a wide range of indications such as: acute and chronic back pain, post-operative support, activity related back pain, disc herniation, stable, non-displaced fractures.
- Removable rigid anterior and posterior panels combined with lateral support panels, provide superior posterior support from T9-S1.



EXOS FORM™ II 631

- The Exos Form II 631 provides a higher degree of spinal support and relief from a wide range of indications
- from acute lower back pain to post-operative support. The Exos Form II 631 provides superior sagittal control and support from T9-S1.
- control and support from T9-S1.



EXOS® FORM II 626

- Provides relief for chronic back pain, post-operative support, activity related back pain, disc herniation.
- Superior sagittal control and support from L-1 to L-5.
- Features removable customisable inserts for additional support.



EXOS® FORM II 621

- Provides stabilisation of the sacroiliac (SI) joint, in case of sacroiliac pain and dysfunction.
- Superior stabilisation of the sacroiliac (SI) joint through direct compression.
- Silicon grip strips prevent the brace from migrating.

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
300621-40	Exos™ Form™ II 621	S/M
300621-60	Exos™ Form™ II 621	L/XL
300626-40	Exos™ Form™ II 626	S/M
300626-60	Exos™ Form™ II 626	L/XL
300631-40	Exos™ Form™ II 631	S/M
300631-60	Exos™ Form™ II 631	L/XL
300637-40	Exos™ Form™ II 637	S/M
300637-60	Exos™ Form™ II 637	L/XL

SIZING INFORMATION - WAIST CIRCUMFERENCE



DESCRIPTION	EXOS FORM™ II 637 AND 631	EXOS FORM™ II 626	EXOS FORM™ II 621
S/M	63cm - 125cm	63cm - 125cm	< 89cm
L/XL	125cm - 152cm	125cm - 152cm	91cm - 140 cm



DONJOY

ISOFORM™ BACK BRACES

Designed to aid superior lower back protection and stabilisation, DonJoy® IsoFORM™ back braces help provide support and pain relief for a comprehensive range of chronic and acute indications.




Designed to aid superior lower back protection and stabilisation, the new DonJoy® IsoFORM™ back braces help provide support and pain relief for a comprehensive range of chronic and acute indications.

ISOFORM FEATURE

Adjustable pull tabs help provide targeted compression and help reduce rotation of brace. The pull tab system independently tightens top and bottom allowing patients to localise compression.

ISOFORM FEATURE

Low profile design is easy to apply and remove. The semi-universal sizing and adjustable belt wings allow clinicians to easily customize the fit for a wide range of patients.





ISOFORM™ LO



INDICATIONS

Low-profile design to help with minor lower back pain

CHARACTERISTICS

- Targeted compression helps support fatigued muscles that spasm and cause pain
- Easy to apply and remove
- Helps provide sagittal control and support from L1 to L5

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1693-2	IsoFORM LO	S/M (Hip 71-122cm)
11-1693-4	IsoFORM LO	L/XL (Hip 122-152cm)
11-1699	Hot/Cold Pack Accessory	-

AVAILABLE ON REQUEST



ISOFORM™ LO+



INDICATIONS

Similar design as the LO option with an anterior panel to help with more trunk support

CHARACTERISTICS

- Acute and chronic low back pain, lower back sprains or strains, osteoporosis, spondylolisthesis

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1694-2	IsoFORM LO+	S/M (Hip 71-122cm)
11-1694-4	IsoFORM LO+	L/XL (Hip 122-152cm)
11-1699	Hot/Cold Pack Accessory	-

AVAILABLE ON REQUEST



ISOFORM™ LSO



INDICATIONS

Trunk stability with motion restriction and pain relief aid

CHARACTERISTICS

- Limited lumbar flexion and extension helps restrict movement
- Targeted compression helps support fatigued muscles that spasm and cause pain
- Low-profile brace for comfortable extended wear
- Easy to apply and remove
- Helps provide superior sagittal control and support from T9-S1

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1695-2	IsoFORM LSO	S/M (Hip 71-122cm)
11-1695-4	IsoFORM LSO	L/XL (Hip 122-152cm)
11-1699	Hot/Cold Pack Accessory	-

AVAILABLE ON REQUEST



ISOFORM™ LSO+



INDICATIONS

Similar design as the LSO option with two lateral extensions to help reduce lateral movement

CHARACTERISTICS

- Helps provide higher degree of spinal support and relief
- Helps provide superior sagittal and coronal control and support from T9-S1

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-1696-2	IsoFORM LSO+	S/M (Hip 71-122cm)
11-1696-4	IsoFORM LSO+	L/XL (Hip 125-152cm)
11-1699	Hot/Cold Pack Accessory	-

AVAILABLE ON REQUEST

DONJOY

ISOFORM TLSO



NEW



INDICATIONS

Suitable for people with compression fractures of the thoracic spine, osteoporosis, thoracic mechanical back pain and kyphosis

CHARACTERISTICS

- Anterior panel - helps restrict flexion and trunk rotation and supports T3-S1
- Posterior panel – rigid paneling helps restrict extension, trunk rotation and supports T3-S1
- Lateral Extensions - Removable semi-rigid supports help reduce lateral movement
- Belt Wings – with removable semi-rigid paneling encompasses the torso for customizable support
- Easy to apply and remove with pivoting quick-release buckles
- Over-the shoulder or underarm strap configuration options

ORDERING INFORMATION

PART NUMBER	HIP CIRCUMFERENCE	SIZE
11-1698-2	71 - 122cm	S/M
11-1698-4	122 - 152cm	L/XL
11-1699	Hot/Cold Pack Accessory	
11-1700	Small Sternal Bar 9"	XS/S
11-1701	Large Sternal Bar 15"	XL/XXL

DONJOY

ISOFORM POSTURAL EXTENSION



NEW



INDICATIONS

Suitable for people with thoracic spine conditions such as kyphosis, osteoporosis, acute pain, mechanical back pain, compression fractures, strain and mild scoliosis.

CHARACTERISTICS

- Adjustable Belt Wings – accommodates different shapes for a custom fit
- Paneling – provides durable rigid support
- Dual Pulley System – provides targeted compression with independent top and bottom compression
- Flexible Pull Tabs – for ease of use with adjustable lace length feature
- Donning Mitt – assists with ease of application
- Semi-Universal – comfortably fits waist sizes from 28-60"
- Hexagonal Spacer Mesh – durable and breathable materials for comfort
- Optional Hot/Cold Pack – for additional relief
- EVA lined Panels – provides comfortable support
- Enclosed Lacing

ORDERING INFORMATION

PART NUMBER	HIP CIRCUMFERENCE	SIZE
11-1697-2	71 - 122cm	S/M
11-1697-4	122 - 152cm	L/XL
11-1699	Hot/Cold Pack Accessory	

Available on request



DONJOY
ACTISTRAP



INDICATIONS

Acute, subacute or chronic lower back pain, lumbar sciatica.

CHARACTERISTICS

- Strong compressive elastic material to enhance posture control with double bilateral straps to adjust compression level
- Four rigid dorsal stays aid stabilisation of the sacro-lumbar area and lordosis correction
- Anatomical design for improved comfort and fit
- Two flexible abdominal stays
- Brace height 26cm (10")

ORDERING INFORMATION

PART NUMBER	WAIST CIRC.	SIZE
82-0231-1	56cm - 68cm (22" - 26.5")	XS
82-0231-2	68cm - 82cm (26.5" - 32")	S
82-0231-3	82cm - 98cm (32" - 38.5")	M
82-0231-4	98cm - 116cm (38.5" - 45.5")	L
82-0231-5	116cm - 136cm (45.5" - 53.5")	XL
82-0231-6	136cm-156cm (53.5" - 61.4")	XXL



DONJOY
**COMFORTFORM
BACK**



INDICATIONS

Moderate compression and support for injuries of the lower back and stabilising support to help prevent spinal overload.

CHARACTERISTICS

- Durable cool mesh elastic for all day wear
- Molded back compression pad aligns and fits into lower lumbar region for maximum support
- Double pull straps help promote proper spinal alignment
- Low profile design for maximum patient comfort and proper fit
- Latex free

ORDERING INFORMATION

PART NUMBER	WAIST CIRC.	SIZE
81-89352	50cm - 63cm (20"-25")	XS
81-89353	63cm - 76cm (25"-30")	S
81-89355	76cm - 86cm (30"-34")	M
81-89357	86cm - 96cm (34"-38")	L
81-89358	96cm - 107 cm (38"-42")	XL
81-89359	107cm - 142cm (42"-56")	XXL
81-89359-10	142cm - 157cm (56"-62")	XXXL



DONJOY
IMMOSTRAP™



INDICATIONS

Acute or subacute lower back pain, lumbar arthrosis, lumbar sciatica, post disc herniation.

CHARACTERISTICS

- Quadruple bilateral straps to adjust compression level
- Independent superior and inferior compression adjustment for optimal anatomical fit
- Ergonomic front panel and a strong compressive elastic material enhances posture control
- Four semi rigid dorsal stays for stabilisation of the sacro-lumbar area
- Rolling buckles for effortless tightening
- Brace height 26cm (10")

ORDERING INFORMATION

PART NUMBER	WAIST CIRC.	SIZE
82-0229-1	56cm - 68cm (22" - 26.5")	XS
82-0229-2	68cm - 82cm (26.5" - 32")	S
82-0229-3	82cm - 98cm (32" - 38.5")	M
82-0229-4	98cm - 116cm (38.5" - 45.5")	L
82-0229-5	116cm - 136cm (45.5" - 53.5")	XL



DONJOY
MYBABYSTRAP



INDICATIONS

- Lower back pain, Pelvic pain, Sacroiliac pain, Postural instabilities
- Discomfort due to stretching and expansion of the abdomen, Back support after birth

CHARACTERISTICS

- Modular design is adaptable to all stages of pregnancy and also after birth. Easy to adjust
- Adjustable back pad for additional support and compression in the back area
- Four posterior rigid stays for extra support in the lumbar area
- Bilateral straps for adjustable relief
- Anatomical shape follows the natural contours of a pregnant body for an improved fit, with no pressure on the abdomen
- Thin, soft and breathable fabric for a comfortable fit
- Ergonomic front panel
- Easy-to-use with simple, hook-loop closure

ORDERING INFORMATION

PART NUMBER	WAIST CIRC.
82-0276	80 - 150cm

Measurement taken during pregnancy

DONJOY

BODYSTRAP



INDICATIONS

Acute, sub-acute or chronic lower back pain. Return to daily living activities, lightweight tasks.

CHARACTERISTICS

- Strong compressive elastic material
- 4 Flexible dorsal stays
- Bilateral elastic straps
- Finger loop for ease of application
- Ergonomic front panel provides optimal comfort while sitting or standing
- Low profile design

ORDERING INFORMATION

PART NUMBER	WAIST CIRCUMFERENCE	SIZE
82-0260-1	56 cm - 68 cm	XS
82-0260-2	68 cm - 82 cm	S
82-0260-3	82 cm - 98 cm	M
82-0260-4	98 cm - 116 cm	L
82-0260-5	116 cm - 136 cm	XL

DONJOY

SPINOSTRAP



INDICATIONS

DonJoy Spinostrap is designed to provide support and stabilisation and partial immobilisation (restriction of movement) of the spine.

CHARACTERISTICS

- Semi-rigid exoskeletal dorsal frame to provide support and aid relief of the thoracic spine, aluminum core can be molded to fit the individual curvature of the back
- Adjustable lumbar belt uses accessible strap tabs and finger pockets to provide adjustable support for the spine while leaving the abdomen free. Fixes easily using hook and loop fastening
- Adjustable shoulder straps use fixable hook-and-loop strap tabs to provide adjustable tension, aiding correction of kyphotic posture by pulling the shoulders back and straightening the back
- PVC and Latex free

ORDERING INFORMATION

PART NUMBER	SPINE LENGTH C7-S1	SIZE
82-0057-1	42cm-52cm	S
82-0057-2	52cm-62cm	M
82-0057-3	62cm-72cm	L

Available on request

NEW



DONJOY

OSTEOSTRAP



INDICATIONS

- Osteoporosis
- Stable vertebral fractures
- Dorsal hyper-kyphosis
- Osteoarthritis
- Traumatic fractures
- Postural dysfunction associated with Parkinson's disease
- Muscular weakness

CHARACTERISTICS

- Reverse shoulder straps allow adjustable tension
- Semi-rigid aluminium dorsal frame can be molded to fit the individual curvature of the back
- Adjustable lumbar belt uses accessible finger pockets and strap tabs


ORDERING INFORMATION

COMPLETE BACK BRACE

PART NUMBER	WAIST CIRC.	SPINE LENGTH	SIZE
82-0049-S1	68 - 82cm	43 - 49cm	S1
82-0049-M1	82 - 98cm	43 - 49cm	M1
82-0049-M2	82 - 98cm	50 - 57cm	M2
82-0049-L2	98 - 116cm	50 - 57cm	L2
82-0049-L3	98 - 116cm	58 - 67cm	L3

SHOULDER & ARM

ULTRASLING PRO



SHOULDER & ARM BRACING

THE BEST FIT GUIDE

SHOULDER		
SMART	STANDARD	VALUE
UltraSling Quadrant	Ultrasling PRO	Sully Shoulder Stabiliser
-	Ultrasling 3	Clinical Shoulder Immobiliser
-	Humeral Fracture Brace/Over Shoulder	Quickfit shoulder immobiliser
-	Humeral Stabilising System (HSS)	-
-	Cromax	-

ARM		
SMART	STANDARD	VALUE
X-Act ROM Elbow	Healthcare Bear Arm Sling	Pneumatic Armband
Epiforce Strong	-	-

COLD THERAPY SHOULDER OPTIONS

For more details see pages 80-89



COLD THERAPY ARM OPTIONS

For more details see pages 80-89



DONJOY

SULLY® SHOULDER STABILISER



INDICATIONS

Shoulder immobilisation and controlled range of motion for anterior, multi-directional, inferior, and posterior instabilities, rotator cuff deceleration, shoulder separations, and muscle strains.

CHARACTERISTICS

- Stabilises, assists and restricts movement of the shoulder post-operatively
- Can pull the humerus posteriorly as it restricts external rotation and abduction
- Elastic straps attach with Velcro at any point, in any direction for improved patient fit
- Smooth, controlled assistance or restraint follows the natural movement of the muscles and joints
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	CHEST CIRC.	BICEP CIRC.	SIZE
11-0525-1	86.5cm - 101.5cm	22.5cm - 26.5cm	XS
11-0525-2	91.5cm - 101.5cm	28cm - 33cm	S
11-0525-3	96.5cm - 106.5cm	33cm - 38cm	M
11-0525-4	106.5cm - 116.5cm	35.5cm - 43cm	L
11-0525-5	116.5cm - 129.5cm	38cm - 43cm	XL

Measurement taken at bicep in non-flex position



PROCARE HUMERAL FRACTURE BRACE/OVER THE SHOULDER



DONJOY HUMERAL STABILISING SYSTEM (HSS)



INDICATIONS

Indicated for management of mid-shaft and distal third humeral fractures.

CHARACTERISTICS

- Pre-formed orthopaedic grade thermoplastic cast-like construction for the humerus
- Hook and loop straps for easy adjustment
- Deltoid cap extension, providing a greater alignment control and prevents slippage
- Supplied with two double thickness stockinettes
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	BICEP CIRC.	LENGTH	SIZE
79-97953	20.5cm-28cm	25.5cm	S
79-97955	25.5cm-33cm	26cm	M
79-97957	30.5cm-38cm	28.5cm	L
79-97958	35.5cm-43cm	32.5cm	XL

INDICATIONS

Immobilisation of stable midshaft humeral fractures, stable supracondylar fractures and biceps/triceps repairs.

CHARACTERISTICS

- Lightweight design offers three point stabilisation
- 900g weight for axial traction
- Full circumference shells and padding with the ability to break down the upper arm cuff to suit all patients
- Alternative to serial casting

ORDERING INFORMATION

PART NUMBER	ARM	SIZE
11-0113-9-13066	Left	U
11-0112-9-13066	Right	U

DONJOY
CROMAX



INDICATIONS

Clavicular fractures or dislocations

CHARACTERISTICS

- Foam/nylon construction with contact closure for easy adjustment
- Movable D-rings allow for improved fit and adjustment

ORDERING INFORMATION

PART NUMBER	SIZE
H-CR122	Adult
H-CR120	Paediatric



PROCARE
**CLINICAL SHOULDER
IMMOBILISER**



INDICATIONS


Preventing shoulder rotation, post trauma after shoulder subluxations and for over use shoulder injuries.

CHARACTERISTICS

- High quality, durable poly cotton sling
- Contact closure adjustment on the shoulder and body straps
- Shoulder pad to provide additional comfort
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	DEPTH x LENGTH	SIZE
79-84012	15cm x 28cm	XS
79-84013	18cm x 33cm	S
79-84015	20.5cm x 38cm	M
79-84017	23cm x 46cm	L
79-84018	25.5cm x 51cm	XL



PROCARE
**HEALTHCARE BEAR
ARM SLING**



INDICATIONS

Injuries to the arm, wrist or hand of paediatric patients.

CHARACTERISTICS

- Made from a cotton/poly envelope with a cotton web strap with quick release buckle for easy application
- Paediatric friendly, bear print design on the sling
- Suitable for use in conjunction with a cast
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	DEPTH x LENGTH	SIZE
79-99132	11.5cm x 28cm	XS / Paediatric
79-99133	16.5cm x 33cm	S / Child

DONJOY

ULTRASLING QUADRANT™



INDICATIONS

Post-operative immobilisation in 15° of abduction or greater for shoulder hemiarthroplasty, total shoulder replacement, shoulder reconstruction, superior labral repair and shoulder debridement, protocols involving stabilisation following fractures to the proximal humerus, elbow and forearm, Bankart procedures, rotator cuff repair, biceps tendon surgery, shoulder dislocation and elbow ligament and tendon procedures.

CHARACTERISTICS

- TempGuard™ - Cool to touch, anti-microbial material and optional padded shoulder strap for extended duration of use
- Optimal versatility provides four angles of immobilisation in one brace; internal, external, abduction or adduction
- ProTherapy™ - Precision based system which allows the patient to flex and extend their arm, promoting therapeutic movement throughout recovery
- Elbow opening to eliminate pressure on the olecranon bursa
- Patient friendly, single handed application for patient compliance and optional shoulder strap
- Optional extension belt available separately

ORDERING INFORMATION

PART NUMBER	ARM	WAIST CIRCUMFERENCE
11-0211	Left	61cm - 114cm
11-0210	Right	61cm - 114cm
25-0212	Extension Belt	114cm - 152cm

Measurement taken at waist circumference 1" (2.5cm) above belly button

DONJOY

ULTRASLING® III



INDICATIONS

Immobilisation for rotator cuff repairs, capsular shifts, Bankart repairs, glenohumeral dislocations/subluxation and soft tissue repairs/strains.

CHARACTERISTICS

- Easy open front panel allowing forearm exercise and exercise ball to stimulate circulation
- Moisture wicking mesh fabric for greater comfort to reduce risk of secondary infections, particularly in longer term use
- Allows arm to be positioned in a variety of positions dependent on post-operative requirements (either more posterior or anterior)
- Quick release buckles make it easy to take on and off for exercise/therapy
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	MEASUREMENT	SIZE
11-0449-2	<28cm	S
11-0449-3	28cm - 33cm	M
11-0449-4	33cm - 38cm	L
11-0449-5	>38cm	XL

Measurement taken from elbow crease to base of index finger

DONJOY

ULTRASLING PRO™



INDICATIONS

Post-operative rotator cuff repairs. Post-operative Bankart procedures. Capsular shifts. Glenohumeral dislocation/subluxation. Soft tissue strains/repairs

CHARACTERISTICS

- Universally designed to fit the right or left arm
- Low profile, adjustable ring easily slides over the non-injured shoulder, helping to evenly distribute weight and take pressure off neck
- Easy, one-handed application
- Post-op step-down: waist belt straps easily adhere directly to sling when pillow is removed
- UltraSling PRO provide 15° abduction without rotation
- UltraSling PRO ER provides 15° abduction with 15° or 30° internal or external rotation
- UltraSling PRO AB provides 45° or 60° abduction without rotation

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-0447-9	UltraSling PRO	Universal
11-0471-9	UltraSling PRO ER/IR 15°	Universal
11-0472-9	UltraSling PRO ER/IR 30°	Universal
11-0473-9	UltraSling PRO AB	Universal

AIRCAST

QUICKFIT SHOULDER IMMOBILISER



INDICATIONS

Support for post-operative rotator cuff repairs, Bankart procedures, capsular shifts, glenohumeral, dislocation/subluxation, soft tissue strains/repairs. The pillow would be considered when the tissues across the top of the glenohumeral joint need to be in a more relaxed position during healing

CHARACTERISTICS

- Constructed of breathable mesh for a cool, comfortable fit
- Release clip design on the waist allows for easy removal without changing the fit adjustment
- Custom fit by folding the sling onto itself, to adjust to the arm length of the patient
- 15° abduction pillow for proper arm placement in internal rotation
- Thumb loop inside the sling to hold the hand in place

ORDERING INFORMATION

PART NUMBER	MIN LENGTH	MAX LENGTH	DEPTH
06AB	26.5 cm	43 cm	18 cm
06A (Sling Only)	26.5 cm	43 cm	18 cm

Measurement taken from elbow crease to base of index finger

X-ACT ROM ELBOW



INDICATIONS

Fixed or controlled range of motion post injury or post operatively, stable fractures, biceps repair, mild to severe soft tissue injuries.

CHARACTERISTICS

- User-friendly range of motion control with flexion settings from -10° to 120° in 10° increments and extension settings from -10° to 90° in 10° increments
- Quick-lock function from 0° to 90° in 15° increments
- Cuffs telescope independently allowing a precise fit on any patient, adjustability from 36cm to 46cm
- Durable lateral aluminium hinge maximises protection and minimises brace weight
- Unique semi-rigid thermoformed and mouldable liners with anti-migration properties

ORDERING INFORMATION

PART NUMBER	ARM	SIZE
11-9122	Left	U
11-9121	Right	U



AIRCAST
PNEUMATIC
ARMBAND



INDICATIONS

Epicondylitis (tennis elbow), medial epicondylitis (golfer's elbow).

CHARACTERISTICS

- Single pre-inflated aircell concentrates compression on the extensor muscle, providing more support, less constriction
- Latex free
- Fits left or right elbow

ORDERING INFORMATION

PART NUMBER	FOREARM CIRC.	SIZE
81-05AB	20cm - 35.5cm	U



DONJOY
EPIFORCE STRONG



INDICATIONS

Epicondyle/ligament and tendon related trauma/injury
Medial or lateral epicondylitis (e.g. golfer's elbow, tennis elbow)
Elbow pain due to overuse, chronic inflammation

CHARACTERISTICS

- Removable pad stimulates the painful tendon insertion on the inner or outer side of the elbow to help relieve pain and regain mobility
- Adjustable strap for precise application of pressure; increase during activity, decrease during rest
- Slap band for easy one-handed application
- Buckle closure for ease in donning and doffing
- Anatomical design for optimal fit
- Low profile to facilitate use under clothing
- Lightweight for comfortable use every day
- Fits left or right arm

ORDERING INFORMATION

PART NUMBER	SIZE
82-1760	Universal



WRIST & HAND

RESPIFORM



WRIST & HAND

THE BEST FIT GUIDE

WRIST & HAND		
SMART	STANDARD	VALUE
Exos Wrist Range	Comfortform wrist	Quickfit Wrist
Respiform	Comfortform wrist with thumb	Quickfit W.T.O
Respiform +	Rhizoform	Universal wrist-o-prene
-	Ergoform	Universal thumb-o-prene
-	Carpaform	Frog/Baseball splint
-	Polyform	Staxx finger splints
-	Digiform	10" lace up
-	Digiform +	-

WRIST & HAND

COLD THERAPY WRIST & HAND OPTIONS

For more details see pages 80-89



CRYO/CUFF™ HAND



EXOS
WRIST BRACE WITH
BOA®



INDICATIONS

Carpal bone injuries such as lunate, pisiform, or triquetral fractures, scapholunate dislocations, triangular fibrocartilage complex tears, radio-carpal ligament injuries, or minimally displaced or stable distal radius or distal styloid fractures.

CHARACTERISTICS

- The openings around the base of the thumb and across the palmar crease enable range of motion and provide greater comfort without loss of stability
- Detachable nylon strap across the web space provides comfortable, adjustable control distally

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
221-21-1111	221-22-1111	16.5cm - 18.5cm	XXS
221-31-1111	221-32-1111	18.5cm - 20cm	XS
221-41-1111	221-42-1111	20cm - 23cm	S
221-51-1111	221-52-1111	23cm - 25.5cm	M
221-61-1111	221-62-1111	25.5cm - 28cm	L
221-71-1111	221-72-1111	>28cm	XL



EXOS
RADIAL GUTTER
FRACTURE BRACE
WITH BOA®



INDICATIONS

Stabilise fractures to the 2nd and/or 3rd metacarpal bones and the respective metacarpophalangeal joints of the hand. Radial gutter portion allows the caregiver to achieve greater than 70° of flexion at the metacarpal joints of the 2nd and 3rd digits.

CHARACTERISTICS

- Full hook-and-loop release ribbon allows for easier application
- Dorsal hook-and-loop secure tab secures the over-lapping (top) portion of brace for better stability
- Hook-and-loop straps across the web space and around the 2nd and 3rd digits can be attached anywhere to the brace surface for better stability
- Trimmable finger section accommodates desired range of motion

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
326-31-1111	326-32-1111	18.5cm - 20cm	XS
326-41-1111	326-42-1111	20cm - 23cm	S
326-51-1111	326-52-1111	23cm - 25.5cm	M
326-61-1111	326-62-1111	25.5cm - 28cm	L
326-71-1111	326-72-1111	>28cm	XL



EXOS
HUMERAL FRACTURE
BRACE WITH BOA®



INDICATIONS

Stabilise fractures to the mid-shaft of the humerus and some injuries to the more proximal humerus.

CHARACTERISTICS

- Brace can be easily trimmed to the appropriate circumference
- Edge tape provided in the package can be easily applied to the trimmed edge of the brace

ORDERING INFORMATION

PART NUMBER	WIDTH	ARM CIRC.	SIZE
330-40-1111	10.5cm	15.5cm - 40.5cm	S
330-50-1111	12.5cm	23cm - 45.5cm	M
330-60-1111	15cm	28cm - 48cm	L



EXOS
SHORT ARM BRACE
WITH BOA® OPEN
THUMB



INDICATIONS

Fractures to the distal radius or distal ulna, and can be employed for acute injuries as well as post-operatively. Support and protection for major injuries to the wrist.

CHARACTERISTICS

- The open thumb design allows for enhanced ease of fit and patient comfort

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
312-31-1111	312-32-1111	18.5cm - 20cm	XS
312-41-1111	312-42-1111	20cm - 23cm	S
312-51-1111	312-52-1111	23cm - 25.5cm	M
312-61-1111	312-62-1111	25.5cm - 28cm	L
312-71-1111	312-72-1111	>28cm	XL



EXOS
THUMB SPICA
FRACTURE BRACE
WITH BOA®



INDICATIONS

Stabilisation to the thumb, first metacarpophalangeal joint, the first carpometacarpal joint, wrist joint, distal radius and distal ulna. Useful for all injuries to the thumb, such as Bennett's or complex Rolando's fractures, suspected or confirmed navicular (scaphoid) fractures, and for other carpal bone injuries.

CHARACTERISTICS

- Bivalve thumb opening accommodates swelling/dressings
- Extended, trimmable length of the thumb for maximum protection
- Full hook-and-loop release ribbon for easier application

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
315-31-1111	315-32-1111	18.5cm - 20cm	XS
315-41-1111	315-42-1111	20cm - 23cm	S
315-51-1111	315-52-1111	23cm - 25.5cm	M
315-61-1111	315-62-1111	25.5cm - 28cm	L
315-71-1111	315-72-1111	>28cm	XL



EXOS
BOXER'S FRACTURE
BRACE WITH BOA®



INDICATIONS

Stabilise fractures to the 4th and/or 5th metacarpal bones, provides stabilisation to the 4th and 5th interphalangeal and wrist joints. May be used for head/neck fracture of the 4th and/or 5th metacarpal.

CHARACTERISTICS

- Dorsal hook-and-loop secure tab secures the overlapping (top) portion of brace for better stability
- Hook-and-loop fastener straps across the web space and around the 4th and 5th digits can be attached anywhere to the brace surface for better stability
- Trimmable finger section accommodates desired range of motion

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
325-31-1111	325-32-1111	18.5cm - 20cm	XS
325-41-1111	325-42-1111	20cm - 23cm	S
325-51-1111	325-52-1111	23cm - 25.5cm	M
325-61-1111	325-62-1111	25.5cm - 28cm	L
325-71-1111	325-72-1111	>28cm	XL



EXOS
PAEDIATRIC SHORT
ARM BRACE WITH
BOA® - OPEN THUMB



INDICATIONS

Support and protection for major injuries to the wrist. Primarily for fractures to the distal radius or distal ulna and can be employed for acute injuries as well as post-operatively.

CHARACTERISTICS

- Open thumb design allows for enhanced ease of fit

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	COLOUR	SIZE
LEFT	RIGHT			
311-11-1111	311-12-1111	14.5cm - 16.5cm	Black	XXXS
311-21-1111	311-22-1111	16.5cm - 18.5cm	Black	XXS
311-11-2293	311-12-2293	14.5cm - 16.5cm	Spider	XXXS
311-21-2293	311-22-2293	16.5cm - 18.5cm	Spider	XXS
311-11-3285	311-12-3285	14.5cm - 16.5cm	Polka Dot	XXXS
311-21-3285	311-22-3285	16.5cm - 18.5cm	Polka Dot	XXS



EXOS
SHORT THUMB
SPICA II



INDICATIONS

Stabilisation to the first metacarpophalangeal joint and may be used pre-operative, post-operative, post-traumatic as well as for degenerative conditions.

CHARACTERISTICS

- Adjustable thumb circumference accommodates a variety of sizes

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
212-31-1111	212-32-1111	18.5cm - 20cm	XS
212-41-1111	212-42-1111	20cm - 23cm	S
212-51-1111	212-52-1111	23cm - 25.5cm	M
212-61-1111	212-62-1111	25.5cm - 28cm	L
212-71-1111	212-72-1111	>28cm	XL



EXOS
EXTENDED SHORT THUMB SPICA



INDICATIONS

Immobilises the base of the thumb and the interphalangeal joint to stabilise injuries to the ulnar collateral ligament (UCL) (Gamekeeper's thumb), fractures or avulsion injuries to the proximal first phalanx or distal first metacarpal, severe sprains of the thumb, or for arthritic conditions involving the metacarpophalangeal joint, or the interphalangeal joint. May also be used to immobilise the thumb joints following surgical repairs and soft tissue injury.

CHARACTERISTICS

- Bivalve thumb opening accommodates swelling or dressings
- Hook-and-loop contact closure

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
211-31-1111	211-32-1111	18.5cm - 20cm	XS
211-41-1111	211-42-1111	20cm - 23cm	S
211-51-1111	211-52-1111	23cm - 25.5cm	M
211-61-1111	211-62-1111	25.5cm - 28cm	L
211-71-1111	211-72-1111	>28cm	XL



EXOS
LONG THUMB SPICA II WITH BOA®



INDICATIONS

Stabilisation to the thumb in adduction, metacarpophalangeal joint, carpometacarpal joint, distal radius, distal ulna and conditions where immobilisation of the wrist in the sagittal and frontal plane is required. May be used pre-operative, post-operative, post-traumatic as well as for carpal tunnel syndrome.

CHARACTERISTICS

- Increased forearm length for better immobilisation and a Velcro release tab allows the brace to open more fully
- Adjustable thumb strap provides better immobilisation

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
232-31-1111	232-32-1111	18.5cm - 20cm	XS
232-41-1111	232-42-1111	20cm - 23cm	S
232-51-1111	232-52-1111	23cm - 25.5cm	M
232-61-1111	232-62-1111	25.5cm - 28cm	L
232-71-1111	232-72-1111	>28cm	XL



EXOS
HAND BASED ULNAR GUTTER BRACE



INDICATIONS

Stabilises fractures to the 4th and/or 5th metacarpals or phalanges. The ulnar gutter portion of the brace allows the caregiver to achieve greater than 70° of flexion at the metacarpal joints of the 4th and 5th digits.

CHARACTERISTICS

- Web straps with hook-and-loop contact closure
- Maintains free range of motion at the wrist
- Entire brace can be trimmed for desired range of motion

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
125-31-1111	125-32-1111	18.5cm - 20cm	XS
125-41-1111	125-42-1111	20cm - 23cm	S
125-51-1111	125-52-1111	23cm - 25.5cm	M
125-61-1111	125-62-1111	25.5cm - 28cm	L
125-71-1111	125-72-1111	>28cm	XL



EXOS
HAND BASED RADIAL GUTTER BRACE



INDICATIONS


Stabilise fractures to the 2nd and/or 3rd metacarpals or phalanges. Radial gutter portion allows the caregiver to achieve greater than 70° of flexion at the metacarpal joints of the 2nd and 3rd digits.

CHARACTERISTICS

- Web straps with hook-and-loop contact closure
- Maintains free range of motion at the wrist
- Entire brace can be trimmed for desired range of motion

ORDERING INFORMATION

PART NUMBER		HAND CIRC.	SIZE
LEFT	RIGHT		
126-31-1111	126-32-1111	18.5cm - 20cm	XS
126-41-1111	126-42-1111	20cm - 23cm	S
126-51-1111	126-52-1111	23cm - 25.5cm	M
126-61-1111	126-62-1111	25.5cm - 28cm	L
126-71-1111	126-72-1111	>28cm	XL



EXOS SHEET MATERIALS



INDICATIONS

Multiple layered thermoformable material in a single construct ideal for custom fabrication for musculoskeletal stabilisation.

CHARACTERISTICS

- 45.5cm x 61cm sheet size
- Length of UBL edging tape is available as an accessory (Size 2.5cm x 60cm) - pack of 10
- Latex free

ORDERING INFORMATION

PART NUMBER	SHEET THICKNESS
800-40-1111	1.0mm
800-50-1111	1.3mm
800-65-1111	1.5mm

PART NUMBER	DESCRIPTION
800-1111	UBL Edging Tape (Pack of 10)



EXOS BOA ADJUSTABLE TENSIONING KIT



INDICATIONS

Universal tensioning kit for use in combination with Exos Sheet Material for creating custom made splints.

CHARACTERISTICS

- Boa reel and lace dynamic closure system allows for expansion to accommodate any initial swelling

ORDERING INFORMATION

PART NUMBER	SIZE
800-ATK	U



EXOS STRAPS



INDICATIONS

Universal straps for use in combination with Exos Sheet Material for creating custom made splints.

CHARACTERISTICS

- Web material straps with hook/loop contact spots at both ends
- Length 12.5cm
- Supplied in a pack of 10

ORDERING INFORMATION

PART NUMBER	LENGTH
800-50-STRAP	12.5cm



EXOS TOP LOADING OVEN



INDICATIONS

Conveniently sized oven used to heat Exos lower and upper extremity products, to gently mould the brace to the patient's limb.

CHARACTERISTICS

- 230v, UK three pin plug
- Electronic thermostat controls temperature range
- Temperature settings (low 82oC, medium 93oC and high 104oC)
- 0 - 30 minutes timer control

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
800-Oven-UK	Exos Oven

DONJOY

RESPIFORM



INDICATIONS

- Wrist sprains and trauma
- Tendinopathy
- Pre- and post-operative immobilisation

CHARACTERISTICS

- Dedicated for either left or right hand use
- Ergonomic, flexible lightweight design
- Anatomically designed, provides tailored individual fit
- Easy and simple application via adjustable straps and "Stop Strap" system
- Two semi-rigid dorsal stays and an aluminium palmar stay for overall wrist immobilisation
- Advanced Ozone™ material, promotes air flow
- Breathable in general providing good skin care and around-the-clock comfort in different heat and humidity conditions

ORDERING INFORMATION

PART NUMBER		WRIST CIRCUMFERENCE	SIZE
LEFT	RIGHT		
82-96150	82-96160	11cm - 14cm	PED/XS
82-96152	82-96162	14cm - 16cm	S
82-96153	82-96163	16cm - 19cm	M
82-96154	82-96164	19cm - 21cm	L
82-96155	82-96165	21cm - 24cm	XL

DONJOY

RESPIFORM PLUS



INDICATIONS

- Wrist sprains and trauma
- Tendinopathy, rhizarthrosis
- Pre- and post-operative immobilisation

CHARACTERISTICS

- Dedicated for either left or right hand use
- Ergonomic, flexible lightweight design
- Anatomically designed, provides tailored individual fit
- Easy and simple application via adjustable straps and "Stop Strap" system
- Two semi-rigid dorsal stays and an aluminium palmar stay for overall wrist immobilisation
- Advanced Ozone™ material, promotes air flow
- Breathable in general providing good skin care and around-the-clock comfort in different heat and humidity conditions
- Modular designed malleable thumb stay for optimal positioning
- Hook and loop system for easy and secure distal thumb fitting

ORDERING INFORMATION

PART NUMBER		WRIST CIRCUMFERENCE	SIZE
LEFT	RIGHT		
82-96170	82-96180	11cm - 14cm	PED/XS
82-96172	82-96182	14cm - 16cm	S
82-96173	82-96183	16cm - 19cm	M
82-96174	82-96184	19cm - 21cm	L
82-96175	82-96185	21cm - 24cm	XL

DONJOY

DIGIFORM



INDICATIONS

- All indications requiring immobilisation of two fingers
- Sprains, subluxations, dislocations, fractures
- Pre- / post-surgery wrist immobilisation

CHARACTERISTICS

- Easy application via 4 adjustable straps and "Stop Strap" system
- Fits either right or left hand
- Lightweight, low profile design
- Prepositional palmar stay provides tailored support
- Aluminium stay is malleable to adapt to healing protocol
- Removable modular components for ease of washing
- Four removable and adjustable straps provide adjustable compression
- Finger tips protective cover

ORDERING INFORMATION

PART NUMBER	WRIST CIRCUMFERENCE	SIZE
82-96200	11cm - 15cm	PED/XS
82-96202	15cm - 19cm	S/M
82-96203	19cm - 23cm	L/XL

DONJOY

DIGIFORM PLUS



INDICATIONS

- All indications requiring immobilisation of three or four fingers
- Sprains, subluxations, dislocations, fractures
- Pre- / post-surgery wrist immobilisation

CHARACTERISTICS

- Easy application via four adjustable straps and "Stop Strap" system
- Fits either right or left hand
- Lightweight, low profile design
- Prepositional aluminium palmar stay provides tailored support
- Removable modular components for ease of washing
- Finger tips protective cover
- Splint adjustable for combinations of three fingers, or for all four fingers

ORDERING INFORMATION

PART NUMBER	WRIST CIRCUMFERENCE	SIZE
82-96210	11cm - 14cm	PED/XS
82-96212	14cm - 16cm	S
82-96213	16cm - 19cm	M
82-96214	19cm - 21cm	L
82-96215	21cm - 24cm	XL

DONJOY

CARPAFORM



INDICATIONS

- Carpal Tunnel Syndrome
- Tendinopathy, Teno-synovitis

CHARACTERISTICS

- Ergonomic, ultra-lightweight design
- Supportive exoskeleton design for decompression of carpal nerve
- Supports the wrist in neutral position, can be used at night
- Dedicated for either left or right hand use
- Easy and simple to adjust and fit
- Available in three sizes

ORDERING INFORMATION

82-96222		WRIST CIRCUMFERENCE	SIZE
LEFT	RIGHT		
82-96222	82-96232	13cm - 16.5cm	S
82-96223	82-96233	16.5cm - 19.5cm	M
82-96224	82-96234	19.5cm - 25cm	L

DONJOY

POLYFORM WRIST



INDICATIONS

Immobilisation of the wrist/thumb/fingers as part of pain management Pain relief associated with rheumatoid arthritis.

CHARACTERISTICS

- Supportive and comfortable thumb-immobilisation piece
- Easy application via five adjustable straps and "Stop Strap" system
- Removable straps for thumb and finger support
- Adjustable aluminium dorsal stay for improved comfort
- Conformable hand stay for optimal adjustment
- Lightweight (180g), low profile design
- Soft, breathable, latex-free material
- Removable modular components for ease of washing
- Fits either right or left hand
- Available in five sizes (XS to XL)

ORDERING INFORMATION

PART NUMBER	WRIST CIRC.	SIZE
82-96300	11 – 14cm	XS
82-96302	14 – 16cm	S
82-96303	16 – 19cm	M
82-96304	19 – 21cm	L
82-96305	21 – 24cm	XL



DONJOY
RHIZOFORM™



INDICATIONS

Functional malleable thumb brace indicated for rhizarthrosis of the thumb.

CHARACTERISTICS

- Ultra-lightweight supportive exoskeleton structure
- Malleable structure provides tailored support
- onSteam® breathable microfiber lining promotes comfort
- Easy and simple to both adjust and fit

ORDERING INFORMATION

PART NUMBER		WRIST CIRC.	SIZE
LEFT	RIGHT		
82-96242	82-96292	13.5cm -16.5cm	S/M
82-96243	82-96293	16.5cm -20cm	L/XL



DONJOY
ERGOFORM™



INDICATIONS

All indications requiring immobilisation of the thumb joints. Indicated for sprains of the thumb, subluxations, dislocations, rhizoarthrosis.

CHARACTERISTICS

- Easy application via adjustable strap and “Stop Strap” system
- Malleable aluminium thumb stay
- Fits left or right hand
- Advanced Ozone™ material promotes airflow
- Lightweight, low-profile design

ORDERING INFORMATION

PART NUMBER	WRIST CIRC.	SIZE
82-96190	12cm - 16cm	0
82-96191	14cm - 18cm	1
82-96192	18cm - 22cm	2



DONJOY
COMFORTFORM
WRIST



INDICATIONS

Treatment of sprains, strains, and control of carpal tunnel syndrome symptoms.

CHARACTERISTICS

- Durable, lightweight Lycra lined foam laminate for breathability and patient comfort
- Dual mouldable aluminium stays support the palmar surface
- Contact closure straps

ORDERING INFORMATION

PART NUMBER		WRIST CIRC.	BRACE LENGTH	SIZE
LEFT	RIGHT			
81-97892	81-97882	11.5cm - 14cm	19cm	XS
81-97893	81-97883	14cm - 16.5cm	19cm	S
81-97895	81-97885	16.5cm - 19cm	20cm	M
81-97897	81-97887	19cm - 21.5cm	21.5cm	L
81-97898	81-97888	>21.5 cm	21.5cm	XL



DONJOY
COMFORTFORM™
WRIST WITH
ABDUCTED THUMB



INDICATIONS

Treatment of sprains, strains, scaphoid injuries, carpal tunnel syndrome, Gamekeeper’s thumb and De Quervain’s syndrome symptoms.

CHARACTERISTICS

- Durable, lightweight Lycra lined foam laminate for breathability and patient comfort
- Dual mouldable aluminium stays support the palmar surface and the extensor surface of the thumb
- Contact closure straps

ORDERING INFORMATION

PART NUMBER		WRIST CIRC.	BRACE LENGTH	SIZE
LEFT	RIGHT			
81-87312	81-87302	11.5cm - 14cm	19cm	XS
81-87313	81-87303	14cm - 16.5cm	19cm	S
81-87315	81-87305	16.5cm - 19cm	20cm	M
81-87317	81-87307	19cm - 21.5cm	21.5cm	L
81-87318	81-87308	>21.5 cm	21.5cm	XL



PROCARE
QUICK-FIT WRIST



INDICATIONS

Wrist sprains, strains and Carpal Tunnel Syndrome

CHARACTERISTICS

- Adjustable dorsal stay for controlled range of motion
- Helps prevent wrist flexion, extension or rotation
- Removable palmar stay may be molded for custom fit
- Durable, lightweight material for comfort

ORDERING INFORMATION

PART NUMBER		SIZE
LEFT	RIGHT	
79-87570	79-87560	U
79-87571	79-87561	XL



PROCARE
QUICK-FIT W.T.O



INDICATIONS

Scaphoid injuries, sprains and strains

CHARACTERISTICS

- Comfortable, lightweight nylon/foam laminate construction
- Adjustable dorsal stay to accommodate various sized patients
- Pre-contoured aluminium palmar and thumb stays may be adjusted for proper angle to address wrist and thumb pathologies
- Elastic straps with D-ring closure provide compressive support and comfort

ORDERING INFORMATION

PART NUMBER		SIZE
LEFT	RIGHT	
79-87490	79-87480	U



PROCARE
UNIVERSAL
WRIST-O-PRENE



INDICATIONS

Carpal tunnel syndrome and wrist sprains and strains.

CHARACTERISTICS

- Wraparound neoprene brace provides compression and support
- Removable contoured aluminium palmar stay
- Inside contact closure allows for ease of application and proper fit

ORDERING INFORMATION

PART NUMBER		WRIST CIRC.	SIZE
LEFT	RIGHT		
81-82471	81-82470	<28cm	U



PROCARE
UNIVERSAL
THUMB-O-PRENE™



INDICATIONS

Tendonitis, basal joint arthritis and instability of the carpometacarpal joint.

CHARACTERISTICS

- Neoprene construction provides compressive wrist support with thumb abduction
- Unique design combines two stays, rigid and spiral plus a movable elastic strap to address various pathologies of the thumb
- One size fits most patients
- Fits left or right hand

ORDERING INFORMATION

PART NUMBER	SIZE
79-82700	U



DONJOY
10" LACE UP



INDICATIONS

Ideal for wrist sprains, strains, and carpal tunnel syndrome.

CHARACTERISTICS

- Perforated suede/flannel construction.
- Single pull laced closure for easy application.
- Length: 10" length.

ORDERING INFORMATION

PART NUMBER		WRIST CIRC.	SIZE
LEFT	RIGHT		
79-87232	79-87222	≤14 cm	XS
79-87233	79-87223	14 - 16,5 cm	S
79-87235	79-87225	16,5 - 19 cm	M
79-87237	79-87227	19 - 21,5 cm	L
79-87238	79-87228	21,5 cm +	XL



PROCARE
FROG SPLINTS /
BASEBALL SPLINTS



INDICATIONS

Helps maintain phalangeal joints in position without bulky dressing.

CHARACTERISTICS

- Padded aluminium splint.
- Easily contours for proper fit.
- Foam Padded.

ORDERING INFORMATION

FROG SPLINTS		
PART NUMBER	LENGTH	SIZE
79-71963-1	7 x 6 cm	S
79-71965-1	6.5 x 6.5 cm	M
79-71967-1	9 x 7.5 cm	L

BASEBALL SPLINTS		
PART NUMBER	LENGTH	SIZE
79-72093-1	10 cm	S
79-72095-1	11.4 cm	M
79-72097-1	12.7 cm	L

PROCARE
STAXX PLASTIC
FINGER SPLINTS



INDICATIONS

Repair of ruptured tendons, fingertip and nail bed injuries and mallet finger.

CHARACTERISTICS

- Lightweight plastic splints
- Pack of padded finger straps available separately for securing the splints
- Starter set of clear finger splints available (various sizes)

ORDERING INFORMATION

PART NUMBER	LENGTH	DIAMETER	SIZE
79-72261	5.3cm	1.5cm	1
79-72262	5.2cm	1.7cm	2
79-72263	5.8cm	1.9cm	3
79-72264	6.3cm	1.9cm	4
79-72265	5.9cm	1.9cm	5
79-72266	6.1cm	2.3cm	5.5
79-72267	6.4cm	2.3cm	6
79-72268	6.3cm	2.5cm	7

PART NUMBER	DESCRIPTION
79-72240	Finger Strap
79-72260	Clear Splint Starter Set

HOT/COLD THERAPY



AIRCAST

CRYO/CUFF™ WRAPS & CUFFS

INDICATIONS

The anatomically designed cuffs provide the therapeutic benefits of controlled cold compression to minimise haemarthrosis, swelling and reduce pain.

CHARACTERISTICS

- Cuffs are anatomically designed to cover all major joints of the body
- Cuffs provide maximum coverage for optimum cryotherapy treatment
- Soothing cold covers the affected area reducing swelling and providing pain relief
- Cuffs are simple to apply and remove and may be detached from the cooler at any time for mobility without interruption of treatment



Foot



Ankle



Thigh



Elbow



Shoulder



Back/Hip/Rib



Knee



Hand



Calf

PART NO.	DESCRIPTION	SIZE		
10B01	Foot CC	L		25-43 cm
10A01	Ankle CC	U	Universal	-
11A01	Knee CC	M		45-60 cm
11B01	Knee CC	L	Knee Circ.	50-80 cm
13C01	Calf CC	U	Calf Circ.	36-51 cm
13A01	Thigh CC	L		43-58 cm
13B01	Thigh CC	XL	Thigh Circ.	51-69 cm

PART NO.	DESCRIPTION	SIZE		
16A01	Hand/Wrist CC	U	-	-
15A01	Elbow CC	U	-	-
12A01	Shoulder CC	U	Torso Circ.	81-122 cm
14A01	Back/Hip/Rib CC	U	-	-



DONJOY

ICEMAN® CLEAR3™



INDICATIONS

The IceMan CLEAR³ cold therapy unit helps reduce pain and swelling, speeding up rehabilitation. It provides extended cold therapy for a variety of indications and protocols. It utilises DonJoy's patented semi-closed loop recirculation system, which maintains more consistent and accurate temperatures than other cold therapy units, in a preset configuration.

CHARACTERISTICS

- Sure fit, stackable, locking lid with watertight seal
- Low profile design, small cooler for easy portability and storage
- Transparent cooler design for easy refill indication
- Preset inlet temperature for simple and convenient, plug and play use
- Practical, integrated hose hanger for ease of use
- Compatible with a variety of pads

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
11-1639	IceMan CLEAR ³



DONJOY

ICEMAN® WRAP-ON COLD PAD



INDICATIONS

Reduce pain, swelling and facilitate rehabilitation.

CHARACTERISTICS

- Universal shape, suitable for use on many parts of the body
- Wrap dimensions 27cm x 28.5cm
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	SIZE
11-0680-9-00000	U

AIRCAST

CRYO/CUFF™ GRAVITY COOLER



INDICATIONS

Ideal for post-operative recovery, trauma, athletic training rooms and home use. An integrated approach to cold therapy, combining cold and compression therapy.

CHARACTERISTICS

- Cold water delivered using a simple, gravity fed system
- Cooler holds water and ice needed for 6 to 8 hours of cryotherapy
- Compatible with all Cryo/Cuff cuffs and ArcticFlow wraps

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
2125	Cryo/Cuff Gravity Cooler
REPLACEMENT PARTS	DESCRIPTION
2130	Cryo/Cuff Cooler Tube
2002	Cryo/Cuff Cooler Replacement Insulation Disc

CHATTANOOGA

FLEXIPAC



INDICATIONS

Hot or cold therapy, when applied directly to the affected area, is used to control pain due to inflammation, swelling, multiple sclerosis, muscle strains, sprains and spasms.

CHARACTERISTICS

- Flexible gel pack conforms to any body contour at any temperature
- Microwaveable to provide soothing heat
- Freezable to provide penetrating cold therapy
- 6 month warranty

ORDERING INFORMATION

PART NUMBER	QUANTITY	SIZE
4026	1	3cm x 15cm
4020	1	13cm x 25cm
4029	1	20cm x 35cm

COMPEX



TRAINING

- Improve muscle strength
- Optimise training
- Prevent injuries

PAIN RELIEF

- Reduce muscle tension
- Relieve muscular and joint pains
- Reduce heavy legs
- Treat tendinitis

WELLNESS

- Sculpt your body
- Tone your muscles
- Shape your body
- Relax

RECOVERY

- Reduce muscle soreness
- Recover after training
- Reviving massage



SWISS TECHNOLOGY



* Scientific studies available on compex.com



Equipped with unique MI (muscle intelligence) technology, Compex® stimulators adapt to each of your muscles to provide the most effective electrostimulation possible, while offering greater comfort and better performance.

MI-SCAN Scans the muscle and automatically adjusts stimulator parameters to your physiology

MI-TENS Makes it easier to adjust stimulation levels for pain relief programs to ensure optimal results

MI-AUTORANGE A brand new feature that enables your Compex module to determine the most appropriate level of stimulation automatically

MI-ACTION Optimises session effectiveness by combining voluntary contractions with muscle stimulated contractions, giving you total control over your workout

MI-RANGE Indicates the optimal level of stimulation to use for Recovery and Massage programs



COMPEX
SP 8.0

INDICATIONS

Our top of the range wireless stimulator designed for athletes who train and compete on a regular basis. SP 8.0 includes our brand new MI-Autorange feature that continually brings you better results

CHARACTERISTICS

- Recover faster
- Enhance performance, strength and endurance
- Helps avoid injury and manage pain
- 40 programs

ORDERING INFORMATION

PART NUMBER

2539116



COMPEX
SP 6.0

INDICATIONS

For dedicated athletes, SP 6.0 includes all the electrostim functions, but with wireless features for freedom of movement. Its easy-to-use color screen includes electrode placement

CHARACTERISTICS

- Improves strength and restores muscle volume
- Helps relax and recover faster
- 30 programs

ORDERING INFORMATION

PART NUMBER

2538116



COMPEX
SP 4.0

INDICATIONS

Includes all the basic electrostimulation features and MI-technology, SP 4.0 also offers rehabilitation programs that can be helpful for athletes with frequent injuries who practice their activity 3 times per week. WIRED TECHNOLOGY


CHARACTERISTICS

- Improves strength and restores muscle volume
- Helps relax and recover faster
- 30 programs

ORDERING INFORMATION

PART NUMBER

2536116



COMPEX
MINI SET

INDICATIONS

Fitness, recovery, and pain relief.

CHARACTERISTICS

- Wireless, mobile-connected muscle stim device. Connecting to your smartphone via Bluetooth, the Compex Mini wireless muscle stim device and TENS machine was created to help everybody embrace their fitness, recovery and pain management, whether on the go, in the gym or at home
- 6 programs

ORDERING INFORMATION

PART NUMBER

INTL-CX192WI04

COMPEX

FIXX 2.0 MASSAGER



INDICATIONS

The Compex Fixx™ 2.0 massager helps warm up muscles before a workout and relieves stiffness after a hard session. Five powerful intensity settings, paired with a high stall force, give the intense pressure needed to relieve the toughest spots safely and comfortably. Interchangeable tips (included) and a rotating head help tailor you massage to the specific needs of different muscle groups. Removable and rechargeable batteries, a compact design, and quiet motor help provide relief anywhere and everywhere.

CHARACTERISTICS

- 5 interchangeable tips
- LED screen display with power and battery level
- 3 arm position options
- 5 speeds: 16, 25, 33, 41 and 50 strokes per second
- Amplitude: 16 mm
- Weight (with battery and tip): 1.3 kg
- Battery life: 2 hours
- Quiet motor: Less than 60 dB
- Powerful motor: Twice the power of Fixx 1.0
- Includes: 1 Compex Fixx 2.0 massager · 1 premium travel case
- 1 lithium ion battery · 1 power adapter · 1 charging dock · 5 tips

ORDERING INFORMATION

PART NUMBER

INTL-CX202WL01

COMPEX

AYRE



INDICATIONS

Sore or aching muscles, pain and swelling of the muscles. Best for aiding in the relief of sore or aching muscles from exercise or recreational activities, helping you to be ready to go at full speed the next day.

CHARACTERISTICS

- 2 modes: Inflate the four chambers one by one or all four chambers at the same time
- Powerful compression: Up to 120 mmHg
- Battery life: Up to 3 hours and it may be used while it is charging
- Ultra portable: No tethers to a power outlet during compression
- Includes: 2 compression boots · 1 premium travel case · 1 charger

ORDERING INFORMATION

PART NUMBER	HEIGHT	SIZE
INTL-CX201FB01	<183 cm	S/M
INTL-CX201FB-L-XL	>184 cm	L/XL



VASCULAR

AIRCAST

VENAFLOW® ELITE



INDICATIONS

VenaFlow is the only VTE prophylaxis compression device that is proven to mimic ambulation, thus reducing VTE rates by 50% vs other devices. This is the best technology to protect patients and healthcare practices from the rising risk of the VTE crisis

- Prophylaxis for deep vein thrombosis
- VenaFlow Elite System combines normal inflation and graduated sequential compression to accelerate venous velocity and enhance fibrinolysis
- For the superior emptying of veins, VenaFlow Elite also incorporates asymmetric compression

CHARACTERISTICS

- Low profile, lightweight design
- Compliance monitor/alarm
- Telescoping bed hanger for custom adjustability
- Automatic Cuff Detection - system identifies and adjusts pressure based on cuffs attached
- Single Leg Operation Button - one push of a button selects single leg operation
- Preset pressures and alarms - no adjustments necessary
- One pump for calf, thigh and foot cuffs provides for ease of use and minimises inventory
- Breathable and comfortable calf, foot and thigh cuffs available

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
30BI-S	VenaFlow Elite, S Mode



AIRCAST
VENAFLOW® ELITE
THIGH CUFF

PART NUMBER	QUANTITY	SIZE
3045	Pair	U



AIRCAST
VENAFLOW® ELITE
FOOT CUFF

PART NUMBER	QUANTITY	SIZE
3046	Pair	U



AIRCAST
VENAFLOW® ELITE
CALF CUFF

PART NUMBER	QUANTITY	MAX CALF CIRC.	SIZE
3040	Pair	48cm	U
3042	Pair	56cm	XL
3043	Pair	76cm	XXL (Bariatric)



AIRCAST
VENAFLOW® TUBING

PART NUMBER	TUBE LENGTH
3008XL	2.6m (8.5')

AIRCAST
VENAPURE®



VenaPure anti-embolism stockings have been designed with Venous Thromboembolism (VTE) prevention in mind. Enovis has been able to use the years of clinical experience gained in the arena of VTE prophylaxis to combine all of these positive attributes into VenaPure.

INDICATIONS

For use in the prophylaxis of deep vein thrombosis and pulmonary embolism in hospitalised patients and for use in other conditions as directed by the clinician.



CHARACTERISTICS

- Graduated compression according to Sigel profile (18mmHg at the ankle, decreasing to 8mmHg at the thigh) to minimise venous stasis and venous dilation
- Inspection toe opening on top of the foot to minimise risk of pressure point and ensure that skin integrity is maintained
- Popliteal break is easily identified by stitch change to ensure that pressure is relieved around the popliteal vein
- Pressure relieving panel at the inner thigh to ensure reduced compression of the femoral vein
- Continuous silicone band at upper thigh provides excellent fixation to the skin, without compromising compression
- Inlay circumferential knit maintains correct positioning of stocking every time
- Colour coded toe and top ensure that the correct size is selected according to patient measurement
- Tape measure included in the dispenser to make measuring of the patient easier
- Available in two styles with lengths and sizes to fit all patients

ORDERING INFORMATION KNEE LENGTH

CALF CIRC.	LEG LENGTH	CODE	SIZE
<30.5cm	<41cm	400R-2	S
	>41cm	400L-2	
30.5cm - 38cm	<43cm	400R-3	M
	>43cm	400L-3	
38cm - 44.5cm	<46cm	400R-4	L
	>46cm	400L-4	
44.5cm - 51cm	<46cm	400R-5	XL
	>46cm	400L-5	
51cm - 58.4cm	<46cm	400R-6	XXL
	>46cm	400L-6	
58.5cm - 66cm	<46cm	400R-7	XXXL
	>46cm	400L-7	

ORDERING INFORMATION THIGH LENGTH

THIGH CIRC.	CALF CIRC.	LEG	CODE	SIZE
<30.5cm		<74cm	401S-2	S
		74cm - 84cm	401R-2	
		>84cm	401L-2	
<63.5cm	30.5cm - 38cm	<74cm	401S-3	M
		74cm - 84cm	401R-3	
		>84cm	401L-3	
	38cm - 44.5cm	<74cm	401S-4	L
		74cm - 84cm	401R-4	
		>84cm	401L-4	
63.5cm - 81.5cm	44.5cm - 55cm	<74cm	401S-5	XL
		74cm - 84cm	401R-5	
		>84cm	401L-5	
81cm - 91cm	55cm - 66cm	<74cm	401S-6	XXL
		74cm - 84cm	401R-6	
		>84cm	401L-6	





DR. COMFORT
ANNIE BLACK



INDICATIONS

Bunions, Charcot, Diabetes, Oedema, Hammertoes,
High instep, Internal brace.

ORDERING INFORMATION

PART NUMBER	SIZE
4510	6-11 (M & W)



DR. COMFORT
LANA BLACK



INDICATIONS

Bunions, Plantar fasciitis

ORDERING INFORMATION

PART NUMBER	SIZE
110910	5-11 (Combined MW)



DR. COMFORT
FLUTE BLACK



INDICATIONS

Diabetes

ORDERING INFORMATION

PART NUMBER	SIZE
2315	6-12 (M & W)



DR. COMFORT
DIANE BLACK



INDICATIONS

Diabetes, hammertoe, plantar fasciitis, stability

ORDERING INFORMATION

PART NUMBER	SIZE
10610	6-12 (M & W)



DR. COMFORT
MALORY BLACK



INDICATIONS

Diabetes, Hammertoe

ORDERING INFORMATION

PART NUMBER	SIZE
10210	6-12 (M & W)

SIZING CHART

WOMEN					
EU	UK	US	EU	UK	US
34	1.5	4	38	5.5	8
34.5	2	4.5	38.5	6	8.5
35	2.5	5	39	6.5	9
35.5	3	5.5	39.5	7	9.5
36	3.5	6	40	7.5	10
36.5	4	6.5	40.5	8	10.5
37	4.5	7	41	8.5	11
37.5	5	7.5	42	9.5	12

At Dr Comfort we understand the importance of foot care and are committed to helping people get back on their feet. With over 100 years combined experience in foot care, footwear knowledge and podiatric training, our team of educators know a thing or two about proper foot health.

Did you know? Revolutionizing the way you walk, our Stability outsole is designed to help provide enhanced support and prevent overpronation. Athletically-inspired shoes that are bold, comfortable and performance-focused. (PDAC approved).

- Wider sole helps provide enhanced stability
- Increased width of toe box provides extra comfort and accommodation
- EVA and rubber outsole provides a layer under foot for added support and comfort
- Firm medial rearfoot helps prevent overpronation



DR. COMFORT
PROTECTOR BLACK



INDICATIONS

Diabetes, ankle stability. Safety boot has a steel cap in the toe box

ORDERING INFORMATION

PART NUMBER	SIZE
9310	8-12 (M, W & XW)



DR. COMFORT
**PERFORMANCE
BLACK/GREY**



INDICATIONS

Bunions, Diabetes, hammertoe, plantar fasciitis, pronation/supination

ORDERING INFORMATION

PART NUMBER	SIZE
7610	8-12 (M & W)



DR. COMFORT
**PERFORMANCE
BLACK**



INDICATIONS

Bunions, Diabetes, hammertoe, plantar fasciitis, pronation/supination

ORDERING INFORMATION

PART NUMBER	SIZE
7611	8-12 (M & W)



DR. COMFORT
BRIAN BLACK



INDICATIONS

Bunions, Charcot, Diabetes, Oedema, Hammertoes, High instep, Internal brace.

ORDERING INFORMATION

PART NUMBER	SIZE
6410	8-12 (M & W)



DR. COMFORT
CLASSIC BLACK



INDICATIONS

Diabetes

ORDERING INFORMATION

PART NUMBER	SIZE
8410	8-12 (M & W)

SIZING CHART

MEN					
EU	UK	US	EU	UK	US
39	5.5	6	43	9.5	10
39.5	6	6.5	44	10	10.5
40	6.5	7	44.5	10.5	11
40.5	7	7.5	45	11	11.5
41	7.5	8	46	11.5	12
41.5	8	8.5	47.5	12.5	13
42	8.5	9	48	13.5	15
42.5	9	9.5	50	14.5	15



PRODUCT WARRANTY

Our products have a limited warranty from the date of purchase (see table below) against defects in materials and workmanship. Normal wear and tear during the usage of a product is not considered a defect.

PRODUCT FAMILY	DONJOY®	PROCARE®	AIRCAST®
CUSTOM DEFIANCE® STEEL HINGE AND FRAME	Lifetime	N/A	N/A
CUSTOM DEFIANCE® STANDARD HINGE AND FRAME	5 years	N/A	N/A
PREFABRICATED RIGID KNEE BRACE FRAME AND HINGE	1 year	6 months	N/A
ASSOCIATED SOFT GOODS*	6 months	6 months	N/A
ELASTIC	6 months	6 months	6 months
WALKER BOOTS	6 months	6 months	6 months
NEOPRENE/DRYTEX	6 months	6 months	6 months
POST-OP BRACING	6 months	6 months	6 months
COLD THERAPY	6 months	N/A	6 months
ASSOCIATED SOFT GOODS*	6 months	N/A	N/A

* Softgoods accessories for rigid and functional bracing requested beyond the 6 month warranty will be sold separately

TERMS & CONDITIONS

1. PARTIES

- a. The "Company" means Enovis South Africa (Pty) Ltd
- b. The "Customer" means the person, firm or company buying goods from the Company to whom the goods are to be invoiced (unless the context otherwise requires) and any person to whom or for the benefit of whom the goods may be delivered.
- c. The "Goods" means the goods described on our invoice.

2. ORDERS

- a. In the Republic of South Africa, credit account customer orders with a value of less than R 9 500.00 (excluding VAT, net discounts) will incur a delivery charge. All orders received from credit account customers in neighboring countries will be responsible for arranging their own collection of stock from our head office.
- b. All orders are to be forwarded to the email address: orders.zaf@enovis.com
- c. V.A.T. Exemption – Buyers who qualify for exemption from paying V.A.T. should include a valid V.A.T. exemption certificate (or other documentation approved by the Customs and Excise for the purpose) with each order.
- d. Orders to non credit account customers require payment to accompany the order before the Goods can be dispatched.

3. QUOTATIONS

Any estimate or quotation shall comprise an invitation to treat and will remain open for a period of 30 days provided that the Company has not previously withdrawn it.

4. SPECIFICATION

Every effort is made to ensure our catalogue descriptions and illustrations are clear and accurate. All colours and sizes are approximate. The reproduction of colours is as accurate as photographic and printing processes allow. Due to our policy of continual improvement, some products may vary from those described. The company reserves the right to amend the specifications, construction and/or design without prior notice.

5. PRICES

- a. All prices quoted are exclusive of Value Added Tax.
- b. The Company reserves the right to vary the price of the Goods to take account of any increases in the cost of raw materials, manufacture, packing, transport, wages or otherwise.
- c. Prices charged will be effective at the date of receipt of Customer order.

6. PAYMENT

- a. Unless the Company otherwise agrees, settlement is due within 30 days of date of statement to Enovis South Africa (Pty) Ltd at the address below.
- b. The Company reserves the right to charge interest on all overdue payments and outstanding balances at the rate of 3% per month.
- c. The Customer shall not be entitled to withhold payment either wholly or in part once the same has become due except where expressly agreed with the Company.
- d. Non credit account customers and sales to members of the public require either payment of a pro-forma invoice or payment to accompany the order before the Goods will be dispatched.
- e. Any invoice dispute needs to be raised within 20 days from receipt of invoice

7. DELIVERY / CARRIAGE CHARGES

- a. Unless the Company otherwise agrees, delivery shall be made to the delivery point specified in the order.
- b. Stated delivery times are an estimate on the part of the Company and shall not be binding upon the Company which shall not be liable for any costs or damages caused by reason of any delay in delivery.
- c. The Company reserve the right to make delivery by instalments and to tender a separate invoice in respect of each instalment.
- d. Orders received before 3.00pm Monday to Friday will only be dispatched the same day for South African customers with credit accounts that do not have overdue payments. This is subject to products being in stock and does not apply to custom manufactured products and non stock items. Delivery will be on a next day basis with some exceptions. Deliveries to remote areas may take up to 72 hours.
- e. Standard Delivery Charges
 - (i) All deliveries to South African residential addresses will incur a delivery charge.
 - (ii) Timed and Saturday deliveries can be arranged at a cost. This will be detailed at the time of ordering.
 - (iii) Oversized items are subject to additional delivery and handling charges, regardless of the value of the order.
- f. Delivery for custom manufactured products may vary from 1 week to 6 weeks from date of order.
- g. Spare parts, unusual sized products and special order items can take between 15 and 30 working days for delivery.

8. CANCELLATIONS

1. No notice of cancellation of an order shall be valid unless given in writing and agreed to by an authorised representative of the Company.
2. If written notice of cancellation is received by the Company, the Company reserves the right to levy a cancellation charge up to 50% of the amount of the order.
3. If the contract provides for delivery by instalments, each instalment shall be deemed to be the subject of a separate contract, and non-delivery or delay in the delivery of any instalment shall not entitle the Customer to repudiate the contract as regards the balance.

9. PASSING OF TITLE

The legal title to the Goods shall not pass to the Customer until all sums due or payable by the Customer to the Company whether in respect of the Goods or otherwise have been received by the Company.

10. PASSING OF RISK

The risk in the Goods shall pass to the Customer when the Goods are tendered for delivery at the delivery point specified in the order. In the event that the Goods are collected by the Customer from the Company, the risk of the Goods shall pass to the Customer when they are loaded onto the vehicle of the Customer or of his carrier or other agent.

11. WARRANTY

The Company warrants that the Goods shall at the time of delivery correspond to the specifications published by the Company when used for purposes for which goods of that type are ordinarily used. In the absence of any such published specifications applicable to the Goods, the Company warrants that the Goods shall at the time of delivery be of merchantable quality. If the Goods do not conform to that warranty, the Company will, at its option:

- a. Replace the Goods found not to conform to the warranty;
- b. Bring the Goods into conformity with the published specifications of the Company or (as the case may be) merchantable quality; or
- c. Take back the Goods found not to conform to the warranty and refund the appropriate part of the purchase price.

12. FORCE MAJEURE

The Company shall not be liable for any loss or damage caused by circumstances beyond the Company's control, including but not limited to, an act of God, war, civil disturbance, governmental restrictions, import or export regulations, industrial disputes and difficulties in obtaining labour or materials. Should any such event occur, the Company might cancel or suspend this contract without incurring any liability for any loss or damage caused.

13. NOTICE OF DAMAGE, NON-DELIVERY OR WRONG DELIVERY

- a. The Company gives no guarantee of the exact quantity of goods delivered, and in the event that the quantity of Goods delivered does not correspond with the quantity stated on the order the customer shall only be liable to pay for the quantity delivered.
- b. If upon receipt the Goods appear damaged or incomplete, please indicate on the carrier's documentation and notify our Customer Services Department of the situation no later than 3 days after delivery. If in doubt, please sign the carrier's documentation "Unchecked". Should an expected consignment not arrive, please notify our Customer Services Department within 21 days from date of invoice, failing which the company shall be under no liability for non-delivery.
- c. Any surplus goods delivered shall remain the property of the Company, and the Customer shall take all reasonable precautions for the safe custody and protection of such surplus goods until the time of their removal by the Company.
- d. In the event of goods being wrongly delivered, under no circumstances shall any goods be returned to the Company without its prior consent, otherwise liability is not accepted. The Company will issue a Returns Authorisation Number (RAN) and this must be marked on the exterior of every package returned. Returned Goods must be suitably packaged for transportation by the carrier. Goods received without suitable packaging will not be credited. Where the Company has made the error, full credit will be given on uplifted goods in good condition.
- e. The Company is under no obligation to accept goods for return, e.g. in the case of Customer ordering error, surplus stock, etc., but where it agrees to do so the following handling charges will apply:-
 - 30 days of delivery, no restocking fee
 - 31 Days – 365 Days of delivery, a 15% restocking/packaging fee
 - Over 12 months of delivery, no credit will be given.The Customer will arrange carriage and pay all freight charges. The Customer should obtain a Returns Authorisation Number, which should be marked clearly on each package being returned. Only complete units of sale may be returned and all goods are to be in good condition. The Customer will be liable for a 15% restocking fee for any goods returned in a non saleable condition or with missing accessories, regardless of the number of days since delivery, with the balance of the goods being paid for in full by the Customer.
- f. Refused Deliveries – We reserve the right to make a carriage charge plus 15% handling charge on all consignments which are tendered for delivery and refused without justification. Goods returned in a non-saleable condition will be paid for in full by the Customer. Goods should not be returned without the prior agreement of the Company.
- g. Only unused and complete goods may be returned, and they must be undamaged, complete with all accessories and less than 12 months old. Goods not conforming to these criteria will not be accepted for credit. The Company reserves the right to dispose of any non saleable and non creditable returned goods.
- h. Custom manufactured products and non stock products are non returnable.

14. DISPUTES AND SET-OFF

Any liability of the Company under the contract shall be subject to and conditional upon the due performance and observance by the Customer of all its obligations under these conditions. Subject to these conditions, the Customer shall not be entitled to withhold or delay payment or exercise any right or set-off whatsoever and howsoever arising or arisen which might otherwise be available to it.

15. COPYRIGHT

All copyrights shall remain the property of the Company alone unless otherwise expressly agreed in writing with the Customer. All product images remain the property of the Company and may not be reproduced in print or electronic format without the express written permission of the Company.

16. RETAIL/RESELLING

Reselling of Enovis branded products is forbidden without written consent from the Company.

17. PATENT RIGHTS

The Customer warrants that any instructions furnished or given by the Customer shall not be such as will cause the Company to infringe any letters patent or trademark or other intellectual property right in execution of the Customer's order and shall indemnify the Company against all actions, costs, claims and such demands arising therefrom to the extent that such instructions lead to such infringements.

18. GOVERNMENT CHARGES

The Customer will reimburse the Company for any increase in tax or governmental charge, or for any new tax or governmental charge hereafter becoming effective which the Company may be required to pay to any government upon the sale, production, or transportation of the Goods and which has the effect of increasing the cost of the Goods sold hereunder.

19. ASSIGNMENT

This contract is personal to the Customer and may only be assigned by it with the prior written consent of the Company.

20. BANKRUPTCY / LIQUIDATION

If the Customer being an individual (or, when the Customer is a firm, any partner in that firm) shall at any time become bankrupt, or having a receiving order made against them or shall make any composition or arrangement with or for the benefit of the creditors or if the Customer, being a company, an order for winding up is made or a resolution for winding up is passed or if a receiver or manager is appointed on behalf of a creditor or otherwise, the Company shall have the right to stop Goods in transit, to suspend further deliveries and to determine any subsisting contract with the Customer and (in any case where title to the Goods has to be passed to the Customer but payment is due and has not been made) the Company shall have a lien over the Goods and right of resale, provided that any such action shall be without prejudice to any right or claim the Company may have against the Customer.

21. CUSTOMER CARE

In our quest to "Never Stop Getting Better", we may from time to time monitor and record customer communications. This is for internal training purposes only and the confidentiality of our Customers and their transactions will be fully protected at all times.

22. PERSONAL CALLERS

We regret that we are unable to deal directly with personal callers to our office.

23. TERMS AND CONDITIONS

Terms and conditions quoted in this catalogue supersede all previous terms and conditions unless otherwise stated.

24. COMPANY DETAILS

Enovis South Africa (Pty) Ltd
Unit 6, 2 on London, 2 London Circle, Brackengate Business Park, Brackenfell, 7560, South Africa
T: 021 276 2968
Company Registration Number: 2009/010124/07 Vat Number: 4080257027

E-mail: orders.zaf@enovis.com
Facebook: Enovis South Africa
Instagram: @enovis_southafrica



enovis[™]

Unit 6, 2 on London, 2 London Circle, Brackengate Business Park, Brackenfell, Cape Town, 7560, South Africa

enovis.co.za | **021 276 2968**

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

MKT-000-00011973-ZA - Rev E

Copyright © 2024 DJO, LLC, a subsidiary of Enovis Corporation