

enovis™

**DONJOY®**

ROAM™ OA





# DONJOY®

Life's next chapter begins with

## ROAM™ OA

### **THE NEWEST INNOVATION IN UNLOADER BRACES FROM DONJOY®, PART OF ENOVIS™**

Help patients say goodbye to knee pain and hello to ROAM OA™. From picking up their grandkids to teeing off at the driving range, **lightweight, low-profile joint unloading and support** never felt so good—or was so easy to prescribe. ROAM helps **improve mobility** and **provides relief** by unloading the pressure of unicompartmental osteoarthritis or other knee pain.

### **COMFORTABLE, EASY TO USE AND EFFECTIVE**

Ideal for a wide variety of patients who want to **increase activity** and who **enjoy an active lifestyle**. For conservative care or in preparation for surgery, ROAM is easy to use and offers a high level of comfort. Consistent use of an unloader knee brace has been shown to contribute to general daily functioning, reduced pain, and improved mobility, which can lead to **increased activity and greater quality of life**.

### **INTENDED USE/INDICATIONS**

For patients with symptoms, physical findings or radiographic evidence of moderate to severe unicompartmental knee osteoarthritis, post-operative meniscal tears only where unicompartmental unloading is indicated, or varus-valgus knee conditions that may benefit from unicompartmental unloading.



### **CUSTOM-FABRICATED OPTION FOR UNIQUE FIT NEEDS**


A custom option addresses your most challenging fitting needs. Custom-measured with condyle and uprights custom-positioned and optimal height setting for the patient's "best fit," ROAM™ provides **low-profile, maximum comfort** and **mitigates migration** by capturing the calf muscle more effectively. Cuff sizes and strap lengths are tailored to the patient. Extension and flexion stops are pre-built into the brace to fulfill the patient's specific needs.

### **COMFORTABLE, EFFECTIVE AND SLEEK**

- **Ultra-low profile, lightweight, fits under clothing**—all-black styling is modern, discreet and breathable
- **Unloads weight** from knee to **dynamically relieve pain**
- **Dynamic strapping and unique condyle harness** completely **avoids sensitive skin in the popliteal space**, while **providing comforting proprioception** on the opposite side of the knee
- **Flexible cuffs** comfortably cradle and **support thigh and calf**
- Designed for a **wide range of body types** and patient needs
- Soft silicone-enhanced lining **feels great on skin; prevents slippage**—2x Cool Technology fibers leave a lasting, cool feeling on the leg
- **Bendable uprights** on the OTS version allow you to better **match patient anatomy**

### **EVERY BRACE ALLOWS TWO FITMENT OPTIONS**

Uniquely designed bracing that can be fitted for either **left lateral and right medial**, or vice versa fitting **right lateral and left medial**. Resulting in a drastically reduced SKU count. **Less is truly more.**





**PATIENT-FOCUSED; PRESCRIBE FOR USE AND SATISFACTION**

- “Set and forget” straps set by a trained fitter
- Two easy connection points; an audible “snap” of the **secure magnetic buckles** gives patients (including those with eyesight or dexterity challenges) **confidence and assurance**
- Patients can **dial in support and pain relief on demand**; BOA® Fit System makes it easy and clear with a visual indicator



DYNAMIC  
UNLOADING


#### QUALITY & DURABILITY

- **Hand washable** after everyday outings and active adventures
- Backed by **DonJoy®** legacy of **quality, innovation, and clinically proven pain relief**
- Compliance—ROAM is **easy to wear**; put it on, take it off and adjust it in seconds
- Designed for efficiency, ease and a smaller environmental footprint; adjustability means just **a few core sizes** can be fit and customized to **most body shapes and sizes**
- **Only 6 off-the-shelf SKUs** are needed to cover the same range that traditionally requires over 50 SKUs


## FEATURES AND BENEFITS




**EASY MAG-ALIGN CLIPS:** These patented magnetic clips help align the strap connections without any fuss. Just hold the attachment using the directional arrow to slide in place and “click” you’re set.




**SET AND FORGET TECHNOLOGY:** You can don and remove your brace without changing any of your doctor’s settings.



**BOA®\* DIAL FOR JUST THE RIGHT FIT:** Dial up or dial down to set exactly how snug the brace fits.



**SLIM PROFILE HINGE:** The slim hinge design allows the brace to be worn under relaxed clothing for a discreet fit.



**ANTI-MIGRATION TECHNOLOGY:** Soft silicone helps keep the brace in place for steady comfort throughout the day.




**PATENT-PENDING COMFORT STRAPS:** The positioning of the support straps is engineered to help the back of your knee remain free from strap rub.

## ORDERING INFORMATION

PART NUMBER		THIGH CIRC.	KNEE CENTER	CALF CIRC.	SIZE
RIGHT MEDIAL/ LEFT LATERAL	LEFT MEDIAL/ RIGHT LATERAL				
11-1726-1-06000	11-1725-1-06000	33 - 47 cm	30.5 - 35.5 cm	25.5 - 35.5 cm	XS/S
11-1726-3-06000	11-1725-3-06000	47 - 59.5 cm	35.5 - 43.2 cm	35.5 - 45.5 cm	M/L
11-1726-5-06000	11-1725-5-06000	59.5 - 81 cm	43.2 - 58.4 cm	45.5 - 61 cm	XL/XXL/XXXL

PART NUMBER	DESCRIPTION
11-R7222-0-00000	Short Thigh, Short Calf, Medial Right/Lateral Left
11-R7223-0-00000	Short Thigh, Short Calf, Medial Left/Lateral Right
11-R7242-0-00000	Short Thigh, Standard Calf, Medial Right/Lateral Left
11-R7243-0-00000	Short Thigh, Standard Calf, Medial Left/Lateral Right
11-R7422-0-00000	Standard Thigh, Short Calf, Medial Right/Lateral Left
11-R7423-0-00000	Standard Thigh, Short Calf, Medial Left/Lateral Right
11-R7442-0-00000	Standard Thigh, Standard Calf, Medial Right/Lateral Left
11-R7443-0-00000	Standard Thigh, Standard Calf, Medial Left/Lateral Right



Scan code for a complete list of accessories and educational resources.



# DONJOY®

## Asia Pacific

1905, Tower 2, Grand Central Plaza  
Shatin  
HONG KONG  
**T:** 852 31051415  
**E:** marketing.ap@enovis.com

## Australia – New Zealand

PO Box 6057  
Frenchs Forest DC  
NSW 2086  
AUSTRALIA  
**T:** 1300 66 77 30  
**E:** customerservice.au@enovis.com

## Benelux

Kleinhoefstraat 5, bus 39  
Geel 2440  
BELGIUM  
**T:** +32 (0)14 248090  
**E:** orders.bx@enovis.com

## Canada

6485 Kennedy Rd  
Mississauga  
Ontario, L5T2W4  
CANADA  
**T:** 1-866-866-5031  
**E:** canada.orders@enovis.com

## China

A312 SOHO ZhongShan Plaza  
1055 W, West ZhongShan Rd  
Shanghai, 200051  
CHINA  
**T:** (8621) 6031 9989  
**E:** shujun.tang@enovis.com

## France

Centre Européen de Frêt  
3 rue de Bethar  
64990 Mouguerre  
FRANCE  
**T:** +33 (0)5 59 52 86 90  
**E:** sce.cial@enovis.com

## Germany

Bötzingen Straße 90  
79111 Freiburg  
GERMANY  
**T:** +49 (0) 180 1 676 333  
**E:** kundenservice@enovis.com

## Italy

Centro Direzionale Milanofiori  
Strada 1, Palazzo F9  
20090 Assago (MI)  
ITALY  
**T:** +39 02 835 98001  
**E:** vendite@enovis.com

## Nordic

Murmansgatan 126  
212 25 Malmö  
SWEDEN  
**T:** +46 40 39 40 00  
**E:** info.nordic@enovis.com

## South Africa

Unit 6, 2 on London, 2 London Circle  
BrackenGate Business Park  
Brackenfell, 7560  
SOUTH AFRICA  
**T:** +27 (0) 21 276 2968  
**E:** info.southafrica@enovis.com

## Spain

Ctra. D'Esplugues, 225 2ª B  
08940 Cornellà de Llobregat  
Barcelona  
SPAIN  
**T:** +34 93 480 32 02  
**E:** ventas@enovis.com

## UK & Ireland

2000 Cathedral Square  
Cathedral Hill  
Guildford, Surrey  
GU2 7YL  
UK  
**T:** +44 (0)1483 459 659  
**E:** ukorders@enovis.com

## USA

2900 Lake Vista Drive  
Lewisville  
TX 75067  
U.S.A.  
**T:** 760 727 1280  
**E:** customercare@enovis.com

**enovis™**

*Creating Better Together™*

**enovis-medtech.eu**

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

\*BOA® is a registered trademark of Boa Technology, Inc.